




PRODUCT CATALOG 2019

-  *Conical Connection*
-  *Internal Hex*

Edison Medical.





Implants Characteristics	2
--------------------------	---

Conical Connection 4 - 50

Implants - Artic	5
Implants - Fortis	6
Implants - ios3	7
Impression & Lab Components	8
Healing Abutments	11
Temporary Abutments	12
Straight Titanium Abutments	14
Angled Titanium Abutment 15°	15
Angled Titanium Abutment 25°	16
Clicq™ Pression - Titanium Transfer Abutments	17
Pop Clicq™ - Overdenture Attachments	22
Pop Clicq™ - Bar	27
Multi Clicq™ - On 4 / On 6	29
Clicq™ - Multi Bar - On 4 / On 6	37
Zircocam™	39
Zircotec™	41
Casting Abutments / UCLA	42
Digital Solutions CAD/CAM	45
Abutment Gauges	49

Internal Hex 51 - 74

Implants - Lamina	52
Impression & Lab Components	53
Healing Abutments	54
Temporary Abutments	55
Titanium Abutments	56
Clicq™ Pression - Titanium Transfer Abutments	59
Pop Clicq™ - Overdenture Attachments	61
Pop Clicq™ - Bar	63
Ball Attachments	64
Q Clicq™ - Overdenture Attachments	65
Multi Clicq™ - On 4 / On 6	67
Clicq™ - Multi Bar - On 4 / On 6	70
Zircocam™ & Zircotec™	71
Casting Abutments & UCLA	72
Digital Solutions CAD/CAM	73
Abutment Gauges	74

Tools 75 - 86

Implant Drivers	76
Screw Drivers	78
Drills	80
Various Tools	83
Surgical Tool Box - Artic	84
Surgical Tool Box - Fortis & iso3	85
Surgical Tool Box - Lamina	86

Implant's Characteristics

Simplicity!
One implant design
various connections



Build-in Platform Switching

Prevent crestal bone loss, which is fundamental for long-term success and stability. It increases the volume of soft tissue around the implant platform, helping to improve the aesthetic final result. As well as reduces the amount of stress transferred from the abutment to the bone.

Micro Rings

Tapered coronal design and micro rings designed to reduce bone stress and minimize bone resorption to optimize bone and soft tissue volume.



Dual Threads



Tapered implant body & core with dual threads and optimal pitch of 1.0mm for faster implant insertion with better stability.

Initial Stability

Progressive threads with compression properties and drilling blades on the apex enable to adjust the implant position during placement.

Implant Surface

Proper implant surface is the key for optimal Osseointegration.
Our unique process is a combination of Sandblasted Large grit particles and Acid etched.
The roughened micro & nano surface structure improve bone adhesion that significantly increase surface area of the implant body for optimal osseointegration.




Safe Insertion

Curved passive apex allows safe implant insertion without the risk of damage to the surrounding.

TI 6AL 4V-ELI Titanium Alloy for Implants
Made in Germany

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Implant's Characteristics



Platform	-A- Overall height	-B- Thread height	-C- Thread pitch	-D- Thread spacing	-E- Collar height	-F- Tip diameter
$\varnothing 3.1\text{mm}$	10	9.5	8.3	2.2	1.1	2.0
	11.5	11	9.8	2.2	1.1	2.0
	13	12.5	11.3	2.2	1.1	2.0
	15	14.5	13.3	2.2	1.1	2.0
$\varnothing 3.5\text{mm}$ $\varnothing 3.6\text{mm}$	8	8	6.8	2.2	1.1	2.1
	10	9.5	8.3	2.2	1.1	2.1
	11.5	11	9.8	2.2	1.1	2.1
	13	12.5	11.3	2.2	1.1	2.1
	15	14.5	13.3	2.2	1.1	2.1
$\varnothing 4.2\text{mm}$ $\varnothing 4.3\text{mm}$	8	8	6.8	2.2	1.1	2.7
	10	9.5	8.3	2.2	1.1	2.7
	11.5	11	9.8	2.2	1.1	2.7
	13	12.5	11.3	2.2	1.1	2.7
	15	14.5	13.3	2.2	1.1	2.7
$\varnothing 5.0\text{mm}$	8	8	6.8	2.2	1.1	3.2
	10	9.5	8.3	2.2	1.1	3.2
	11.5	11	9.8	2.2	1.1	3.2
	13	12.5	11.3	2.2	1.1	3.2

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Conical Connection



Artic® ACTIVE HEX Conical Connection


NobelActive™ Compatible*

 Narrow Platform

\emptyset	L	3.5mm
8mm		AR-NA-NP3508
10mm		AR-NA-NP3510
11.5mm		AR-NA-NP3511
13mm		AR-NA-NP3513
15mm		AR-NA-NP3515

 Regular Platform

\emptyset	L	4.3mm	5.0mm
8mm		AR-NA-RP4308	AR-NA-RP5008
10mm		AR-NA-RP4310	AR-NA-RP5010
11.5mm		AR-NA-RP4311	AR-NA-RP5011
13mm		AR-NA-RP4313	AR-NA-RP5013
15mm		AR-NA-RP4315	




* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.


STREAMLINE
Conical Connection

Osstem® Hiossen® Compatible*

Regular Platform

	Ø	3.6mm	4.2mm	5.0mm
L	8mm	FR-TS-RP3508	FR-TS-RP4208	FR-TS-RP5008
	10mm	FR-TS-RP3510	FR-TS-RP4210	FR-TS-RP5010
	11.5mm	FR-TS-RP3511	FR-TS-RP4211	FR-TS-RP5011
	13mm	FR-TS-RP3513	FR-TS-RP4213	FR-TS-RP5013
	15mm	FR-TS-RP3515	FR-TS-RP4215	



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



ios[®]3




SPEEDY HEX
Conical Connection

Astra Tech™ TX Compatible*

 Slim Platform

	3.1mm
Ø	
L	
10mm	i3-OS-SP3010
11.5mm	i3-OS-SP3011
13mm	i3-OS-SP3013
15mm	i3-OS-SP3015



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Titanium Implant Analog
Regular Platform

	L	10mm	15mm	12mm	15.5mm	12mm
NobelActive™ Compatible *						IA-NA-RP
Osstem® Hiossen® Compatible *						IA-TS-RP
Dentium® Compatible *						IA-DE-RP
MegaGen® AnyOne Compatible *						IA-AO-RP
MegaGen® AnyRidge Compatible *						IA-AR-RP
Astra Tech™TX Compatible *						IA-OS-WP
Dio Implant® Compatible *						IA-UF-RP
NeoBiotech® Compatible *						IA-IS-RP

Titanium Implant Analog
Narrow Platform

	L	12mm	15mm	
NobelActive™ Compatible *				IA-NA-NP
Osstem® Hiossen® Compatible *				IA-TS-NP
Astra Tech™TX Compatible *				IA-OS-SP


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



 Closed Tray - Slim Body - Ø 3.8mm
Regular Platform

	L	9mm	12mm
<hr/>			
NobelActive™ Compatible *	CT-NA-RPSN	CT-NA-RPLN	
Osstem® Hiossen® Compatible *	CT-TS-RPSN	CT-TS-RPLN	
Dentium® Compatible *	CT-DE-RPSN	CT-DE-RPLN	
MegaGen® AnyOne Compatible *	CT-AO-RPSN	CT-AO-RPLN	
MegaGen® AnyRidge Compatible *	CT-AR-RPSN	CT-AR-RPLN	
Astra Tech™TX Compatible *	CT-OS-WPSN	CT-OS-WPLN	
Dio Implant® Compatible *	CT-UF-RPSN	CT-UF-RPLN	
NeoBiotech® Compatible *	CT-IS-RPSN	CT-IS-RPLN	


 Closed Tray - Ø 4.5mm
Regular Platform

	L	9mm	12mm
<hr/>			
NobelActive™ Compatible *	CT-NA-RPS	CT-NA-RPL	
Osstem® Hiossen® Compatible *	CT-TS-RPS	CT-TS-RPL	
Dentium® Compatible *	CT-DE-RPS	CT-DE-RPL	
MegaGen® AnyOne Compatible *	CT-AO-RPS	CT-AO-RPL	
MegaGen® AnyRidge Compatible *	CT-AR-RPS	CT-AR-RPL	
Astra Tech™TX Compatible *	CT-OS-WPS	CT-OS-WPL	
Dio Implant® Compatible *	CT-UF-RPS	CT-UF-RPL	
NeoBiotech® Compatible *	CT-IS-RPS	CT-IS-RPL	



Included

 Closed Tray - Slim Body - Ø 3.8mm
Narrow Platform

	L	9mm	12mm
<hr/>			
NobelActive™ Compatible *	CT-NA-NPSN	CT-NA-NPLN	
Osstem® Hiossen® Compatible *	CT-TS-NPSN	CT-TS-NPLN	
Astra Tech™TX Compatible *	CT-OS-SPSN	CT-OS-SPLN	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Open Tray Regular Platform					
	L		Included		Included
NobelActive™ Compatible*		OT-NA-RP			
Osstem® Hiossen® Compatible*		OT-TS-RP			
Dentium® Compatible*		OT-DE-RP			
MegaGen® AnyOne Compatible *		OT-AO-RP			
MegaGen® AnyRidge Compatible *			OT-AR-RP		
Astra Tech™TX Compatible*		OT-OS-WP			
Dio Implant® Compatible*		OT-UF-RP			
NeoBiotech® Compatible*		OT-IS-RP			

Open Tray Narrow Platform			
	L		Included
NobelActive™ Compatible*		OT-NA-NP	
Osstem® Hiossen® Compatible*		OT-TS-NP	
Astra Tech™TX Compatible*		OT-OS-SP	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Healing Abutment Titanium

Titanium Healing Abutment Regular Platform

H	3mm 5mm	3mm 6mm	
NobelActive™ Compatible *	HA-NA-RP30N	HA-NA-RP30	
Osstem® Hiossen® Compatible *	HA-TS-RP30N	HA-TS-RP30	
Dentium® Compatible *	HA-DE-RP30N	HA-DE-RP30	
MegaGen® AnyOne Compatible *	HA-AO-RP30N	HA-AO-RP30	
MegaGen® AnyRidge Compatible *	HA-AR-RP30N	HA-AR-RP30	
Astra Tech™TX Compatible *	HA-OS-WP30N	HA-OS-WP30	
Dio Implant® Compatible *	HA-UF-RP30N	HA-UF-RP30	
NeoBiotech® Compatible *	HA-IS-RP30N	HA-IS-RP30	

Titanium Healing Abutment Regular Platform

H	5mm 5mm	5mm 6mm	
NobelActive™ Compatible *	HA-NA-RP50N	HA-NA-RP50	
Osstem® Hiossen® Compatible *	HA-TS-RP50N	HA-TS-RP50	
Dentium® Compatible *	HA-DE-RP50N	HA-DE-RP50	
MegaGen® AnyOne Compatible *	HA-AO-RP50N	HA-AO-RP50	
MegaGen® AnyRidge Compatible *	HA-AR-RP50N	HA-AR-RP50	
Astra Tech™TX Compatible *	HA-OS-WP50N	HA-OS-WP50	
Dio Implant® Compatible *	HA-UF-RP50N	HA-UF-RP50	
NeoBiotech® Compatible *	HA-IS-RP50N	HA-IS-RP50	

Titanium Healing Abutment Narrow Platform

H	3mm 4mm	3mm 5mm	
NobelActive™ Compatible *	HA-NA-NP30N	HA-NA-NP30	
Osstem® Hiossen® Compatible *	HA-TS-NP30N	HA-TS-NP30	
Astra Tech™TX Compatible *	HA-OS-SP30N	HA-OS-SP30	

Titanium Healing Abutment Narrow Platform

H	5mm 4mm	5mm 5mm	
NobelActive™ Compatible *	HA-NA-NP50N	HA-NA-NP50	
Osstem® Hiossen® Compatible *	HA-TS-NP50N	HA-TS-NP50	
Astra Tech™TX Compatible *	HA-OS-SP50N	HA-OS-SP50	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Temporary Abutment Titanium

Temporary Titanium Abutment
 Regular Platform

	12mm	12mm
NobelActive™ Compatible*	TTA-NA-RP	TTA-NA-RPR
Osstem® Hiossen® Compatible*	TTA-TS-RP	TTA-TS-RPR
Dentium® Compatible*	TTA-DE-RP	TTA-DE-RPR
MegaGen® AnyOne Compatible*	TTA-AO-RP	TTA-AO-RPR
MegaGen® AnyRidge Compatible*	TTA-AR-RP	TTA-AR-RPR
Astra Tech™TX Compatible*	TTA-OS-WP	TTA-OS-WPR
Dio Implant® Compatible*	TTA-UF-RP	TTA-UF-RPR
NeoBiotech® Compatible*	TTA-IS-RP	TTA-IS-RPR

Temporary Titanium Abutment
Narrow Platform

	12mm	12mm
NobelActive™ Compatible*	TTA-NA-NP	TTA-NA-NPR
Osstem® Hiossen® Compatible*	TTA-TS-NP	TTA-TS-NPR
Astra Tech™TX Compatible*	TTA-OS-SP	TTA-OS-SPR


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection




Temporary Abutment PEEK

PEEK - Straight Abutment
Regular Platform




NobelActive™ Compatible *	PA-NA-RP15	PA-NA-RP30	PA-NA-RP45
Osstem® Hiossen® Compatible *	PA-TS-RP15	PA-TS-RP30	PA-TS-RP45
Dentium® Compatible *	PA-DE-RP15	PA-DE-RP30	PA-DE-RP45
MegaGen® AnyOne Compatible *	PA-AO-RP15	PA-AO-RP30	PA-AO-RP45
MegaGen® AnyRidge Compatible *	PA-AR-RP15	PA-AR-RP30	PA-AR-RP45
Astra Tech™TX Compatible *	PA-OS-WP15	PA-OS-WP30	PA-OS-WP45
Dio Implant® Compatible *	PA-UF-RP15	PA-UF-RP30	PA-UF-RP45
NeoBiotech® Compatible *	PA-IS-RP15	PA-IS-RP30	PA-IS-RP45

PEEK - Anatomic Angled Abutment -15°
Regular Platform




NobelActive™ Compatible *	PAA-NA-RP1515	PAA-NA-RP1530	PAA-NA-RP1545
Osstem® Hiossen® Compatible *	PAA-TS-RP1515	PAA-TS-RP1530	PAA-TS-RP1545
Dentium® Compatible *	PAA-DE-RP1515	PAA-DE-RP1530	PAA-DE-RP1545
MegaGen® AnyOne Compatible *	PAA-AO-RP1515	PAA-AO-RP1530	PAA-AO-RP1545
MegaGen® AnyRidge Compatible *	PAA-AR-RP1515	PAA-AR-RP1530	PAA-AR-RP1545
Astra Tech™TX Compatible *	PAA-OS-WP1515	PAA-OS-WP1530	PAA-OS-WP1545
Dio Implant® Compatible *	PAA-UF-RP1515	PAA-UF-RP1530	PAA-UF-RP1545
NeoBiotech® Compatible *	PAA-IS-RP1515	PAA-IS-RP1530	PAA-IS-RP1545

PEEK - Straight Abutment
Narrow Platform



NobelActive™ Compatible *	PA-NA-NP15	PA-NA-NP30	PA-NA-NP45
Osstem® Hiossen® Compatible *	PA-TS-NP15	PA-TS-NP30	PA-TS-NP45
Astra Tech™TX Compatible *	PA-OS-SP15	PA-OS-SP30	PA-OS-SP45

PEEK - Anatomic Angled Abutment -15°
Narrow Platform



NobelActive™ Compatible *	PAA-NA-NP1515	PAA-NA-NP1530	PAA-NA-NP1545
Osstem® Hiossen® Compatible *	PAA-TS-NP1515	PAA-TS-NP1530	PAA-TS-NP1545
Astra Tech™TX Compatible *	PAA-OS-SP1515	PAA-OS-SP1530	PAA-OS-SP1545

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Straight Abutment Titanium

Straight Titanium Abutment W/Shoulder Regular Platform

		1.5mm	3mm	4.5mm
H				
NobelActive™ Compatible *	TA-NA-RP15	TA-NA-RP30	TA-NA-RP45	
Osstem® Hiossen® Compatible *	TA-TS-RP15	TA-TS-RP30	TA-TS-RP45	
Dentium® Compatible *	TA-DE-RP15	TA-DE-RP30	TA-DE-RP45	
MegaGen® AnyOne Compatible *	TA-AO-RP15	TA-AO-RP30	TA-AO-RP45	
MegaGen® AnyRidge Compatible *	TA-AR-RP15	TA-AR-RP30	TA-AR-RP45	
Astra Tech™TX Compatible *	TA-OS-WP15	TA-OS-WP30	TA-OS-WP45	
Dio Implant® Compatible *	TA-UF-RP15	TA-UF-RP30	TA-UF-RP45	
NeoBiotech® Compatible *	TA-IS-RP15	TA-IS-RP30	TA-IS-RP45	

Anatomic Straight Titanium Abutment Regular Platform

		1.5mm 9.5mm	3mm 11mm	4.5mm 12.5mm
H				
L				
NobelActive™ Compatible *	TAA-NA-RP15	TAA-NA-RP30	TAA-NA-RP45	
Osstem® Hiossen® Compatible *	TAA-TS-RP15	TAA-TS-RP30	TAA-TS-RP45	
Dentium® Compatible *	TAA-DE-RP15	TAA-DE-RP30	TAA-DE-RP45	
MegaGen® AnyOne Compatible *	TAA-AO-RP15	TAA-AO-RP30	TAA-AO-RP45	
MegaGen® AnyRidge Compatible *	TAA-AR-RP15	TAA-AR-RP30	TAA-AR-RP45	
Astra Tech™TX Compatible *	TAA-OS-WP15	TAA-OS-WP30	TAA-OS-WP45	
Dio Implant® Compatible *	TAA-UF-RP15	TAA-UF-RP30	TAA-UF-RP45	
NeoBiotech® Compatible *	TAA-IS-RP15	TAA-IS-RP30	TAA-IS-RP45	

Straight Titanium Abutment W/Shoulder Narrow Platform

		1.5mm	3mm	4.5mm
H				
NobelActive™ Compatible *	TA-NA-NP15	TA-NA-NP30	TA-NA-NP45	
Osstem® Hiossen® Compatible *	TA-TS-NP15	TA-TS-NP30	TA-TS-NP45	
Astra Tech™TX Compatible *	TA-OS-SP15	TA-OS-SP30	TA-OS-SP45	

Anatomic Straight Titanium Abutment Narrow Platform

		1.5mm 9.5mm	3mm 11mm	4.5mm 12.5mm
H				
L				
NobelActive™ Compatible *	TAA-NA-NP15	TAA-NA-NP30	TAA-NA-NP45	
Osstem® Hiossen® Compatible *	TAA-TS-NP15	TAA-TS-NP30	TAA-TS-NP45	
Astra Tech™TX Compatible *	TAA-OS-SP15	TAA-OS-SP30	TAA-OS-SP45	


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection




Angled Abutment 15° Titanium

Standard Angled Abutment 15° Regular Platform




NobelActive™ Compatible*	TA-NA-RP1500
Osstem® Hiossen® Compatible*	TA-TS-RP1500
Dentium® Compatible*	TA-DE-RP1500
MegaGen® AnyOne Compatible*	TA-AO-RP1500
MegaGen® AnyRidge Compatible*	TA-AR-RP1500
Astra Tech™TX Compatible*	TA-OS-WP1500
Dio Implant® Compatible*	TA-UF-RP1500
NeoBiotech® Compatible*	TA-IS-RP1500

Anatomic Angled Abutment 15° Regular Platform




NobelActive™ Compatible*	TAA-NA-RP1515	TAA-NA-RP1530	TAA-NA-RP1545
Osstem® Hiossen® Compatible*	TAA-TS-RP1515	TAA-TS-RP1530	TAA-TS-RP1545
Dentium® Compatible*	TAA-DE-RP1515	TAA-DE-RP1530	TAA-DE-RP1545
MegaGen® AnyOne Compatible*	TAA-AO-RP1515	TAA-AO-RP1530	TAA-AO-RP1545
MegaGen® AnyRidge Compatible*	TAA-AR-RP1515	TAA-AR-RP1530	TAA-AR-RP1545
Astra Tech™TX Compatible*	TAA-OS-WP1515	TAA-OS-WP1530	TAA-OS-WP1545
Dio Implant® Compatible*	TAA-UF-RP1515	TAA-UF-RP1530	TAA-UF-RP1545
NeoBiotech® Compatible*	TAA-IS-RP1515	TAA-IS-RP1530	TAA-IS-RP1545

Standard Angled Abutment 15° Narrow Platform



NobelActive™ Compatible*	TA-NA-NP1500
Osstem® Hiossen® Compatible*	TA-TS-NP1500
Astra Tech™TX Compatible*	TA-OS-SP1500

Anatomic Angled Abutment 15° Narrow Platform



NobelActive™ Compatible*	TAA-NA-NP1515	TAA-NA-NP1530	TAA-NA-NP1545
Osstem® Hiossen® Compatible*	TAA-TS-NP1515	TAA-TS-NP1530	TAA-TS-NP1545
Astra Tech™TX Compatible*	TAA-OS-SP1515	TAA-OS-SP1530	TAA-OS-SP1545


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection




Angled Abutment 25° Titanium

Standard Angled Abutment 25° Regular Platform




NobelActive™ Compatible*	TA-NA-RP2500
Osstem® Hiossen® Compatible*	TA-TS-RP2500
Dentium® Compatible*	TA-DE-RP2500
MegaGen® AnyOne Compatible*	TA-AO-RP2500
MegaGen® AnyRidge Compatible*	TA-AR-RP2500
Astra Tech™TX Compatible*	TA-OS-WP2500
Dio Implant® Compatible*	TA-UF-RP2500
NeoBiotech® Compatible*	TA-IS-RP2500

Anatomic Angled Abutment 25° Regular Platform




NobelActive™ Compatible*	TAA-NA-RP2515	TAA-NA-RP2530	TAA-NA-RP2545
Osstem® Hiossen® Compatible*	TAA-TS-RP2515	TAA-TS-RP2530	TAA-TS-RP2545
Dentium® Compatible*	TAA-DE-RP2515	TAA-DE-RP2530	TAA-DE-RP2545
MegaGen® AnyOne Compatible*	TAA-AO-RP2515	TAA-AO-RP2530	TAA-AO-RP2545
MegaGen® AnyRidge Compatible*	TAA-AR-RP2515	TAA-AR-RP2530	TAA-AR-RP2545
Astra Tech™TX Compatible*	TAA-OS-WP2515	TAA-OS-WP2530	TAA-OS-WP2545
Dio Implant® Compatible*	TAA-UF-RP2515	TAA-UF-RP2530	TAA-UF-RP2545
NeoBiotech® Compatible*	TAA-IS-RP2515	TAA-IS-RP2530	TAA-IS-RP2545

Standard Angled Abutment 25° Narrow Platform




NobelActive™ Compatible*	TA-NA-NP2500
Osstem® Hiossen® Compatible*	TA-TS-NP2500
Astra Tech™TX Compatible*	TA-OS-SP2500

Anatomic Angled Abutment 25° Narrow Platform



NobelActive™ Compatible*	TAA-NA-NP2515	TAA-NA-NP2530	TAA-NA-NP2545
Osstem® Hiossen® Compatible*	TAA-TS-NP2515	TAA-TS-NP2530	TAA-TS-NP2545
Astra Tech™TX Compatible*	TAA-OS-SP2515	TAA-OS-SP2530	TAA-OS-SP2545

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Professional Kit Regular Platform

NobelActive™ Compatible *	TAP-NA-RPPK
Astra Tech™TX Compatible *	TAP-OS-WPPK
Dio Implant® Compatible *	TAP-UF-RPPK
NeoBiotech® Compatible *	TAP-IS-RPPK

Professional Kit Narrow Platform

NobelActive™ Compatible *	TAP-NA-NPPK
Osstem® Hiossen® Compatible *	TAP-TS-NPPK
Astra Tech™TX Compatible *	TAP-OS-SPPK

Clicq™ Pression Components

Abutment Analog



Medical Peek - Classix

Healing Cap PEEK




Waxing Sleeve



Plastic Transfer



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Starter Kit
Regular Platform

NobelActive™ Compatible *	TAP-NA-RPSK
Astra Tech™TX Compatible *	TAP-OS-WPSK
Dio Implant® Compatible *	TAP-UF-RPSK
NeoBiotech® Compatible *	TAP-IS-RPSK

Osstem® Hiossen® Compatible *	TAP-TS-RPSK
Dentium® Compatible *	TAP-DE-RPSK
MegaGen® AnyOne Compatible *	TAP-AO-RPSK

Starter Kit
Narrow Platform

NobelActive™ Compatible *	TAP-NA-NPSK
Osstem® Hiossen® Compatible *	TAP-TS-NPSK
Astra Tech™TX Compatible *	TAP-OS-SPSK



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Transfer Abutment Standard Kit - Titanium



Straight Transfer Abutment - Short
Regular Platform




NobelActive™ Compatible *	TAP-NA-RP15S	TAP-NA-RP30S	TAP-NA-RP45S
Osstem® Hiossen® Compatible *	TAP-TS-RP15S	TAP-TS-RP30S	TAP-TS-RP45S
Dentium® Compatible *	TAP-DE-RP15S	TAP-DE-RP30S	TAP-DE-RP45S
MegaGen® AnyOne Compatible *	TAP-AO-RP15S	TAP-AO-RP30S	TAP-AO-RP45S
MegaGen® AnyRidge Compatible *	TAP-AR-RP15S	TAP-AR-RP30S	TAP-AR-RP45S
Astra Tech™TX Compatible *	TAP-OS-WP15S	TAP-OS-WP30S	TAP-OS-WP45S
Dio Implant® Compatible *	TAP-UF-RP15S	TAP-UF-RP30S	TAP-UF-RP45S
NeoBiotech® Compatible *	TAP-IS-RP15S	TAP-IS-RP30S	TAP-IS-RP45S

Contents:

Straight / Angled Abutment
Titanium Screw
Plastic Transfer



Straight Titanium Transfer Abutment - Short
Narrow Platform



NobelActive™ Compatible *	TAP-NA-NP15S	TAP-NA-NP30S	TAP-NA-NP45S
Osstem® Hiossen® Compatible *	TAP-TS-NP15S	TAP-TS-NP30S	TAP-TS-NP45S
Astra Tech™TX Compatible *	TAP-OS-SP15S	TAP-OS-SP30S	TAP-OS-SP45S


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection




Transfer Abutment Standard Kit - Titanium

Straight Transfer Abutment - Medium
Regular Platform




NobelActive™ Compatible *	TAP-NA-RP15M	TAP-NA-RP30M	TAP-NA-RP45M
Osstem® Hiossen® Compatible *	TAP-TS-RP15M	TAP-TS-RP30M	TAP-TS-RP45M
Dentium® Compatible *	TAP-DE-RP15M	TAP-DE-RP30M	TAP-DE-RP45M
MegaGen® AnyOne Compatible *	TAP-AO-RP15M	TAP-AO-RP30M	TAP-AO-RP45M
MegaGen® AnyRidge Compatible *	TAP-AR-RP15M	TAP-AR-RP30M	TAP-AR-RP45M
Astra Tech™TX Compatible *	TAP-OS-WP15M	TAP-OS-WP30M	TAP-OS-WP45M
Dio Implant® Compatible *	TAP-UF-RP15M	TAP-UF-RP30M	TAP-UF-RP45M
NeoBiotech® Compatible *	TAP-IS-RP15M	TAP-IS-RP30M	TAP-IS-RP45M

Straight Transfer Abutment - Long
Regular Platform




NobelActive™ Compatible *	TAP-NA-RP15L	TAP-NA-RP30L	TAP-NA-RP45L
Osstem® Hiossen® Compatible *	TAP-TS-RP15L	TAP-TS-RP30L	TAP-TS-RP45L
Dentium® Compatible *	TAP-DE-RP15L	TAP-DE-RP30L	TAP-DE-RP45L
MegaGen® AnyOne Compatible *	TAP-AO-RP15L	TAP-AO-RP30L	TAP-AO-RP45L
MegaGen® AnyRidge Compatible *	TAP-AR-RP15L	TAP-AR-RP30L	TAP-AR-RP45L
Astra Tech™TX Compatible *	TAP-OS-WP15L	TAP-OS-WP30L	TAP-OS-WP45L
Dio Implant® Compatible *	TAP-UF-RP15L	TAP-UF-RP30L	TAP-UF-RP45L
NeoBiotech® Compatible *	TAP-IS-RP15L	TAP-IS-RP30L	TAP-IS-RP45L

Straight Transfer Abutment - Medium
Narrow Platform



NobelActive™ Compatible *	TAP-NA-NP15M	TAP-NA-NP30M	TAP-NA-NP45M
Osstem® Hiossen® Compatible *	TAP-TS-NP15M	TAP-TS-NP30M	TAP-TS-NP45M
Astra Tech™TX Compatible *	TAP-OS-SP15M	TAP-OS-SP30M	TAP-OS-SP45M

Straight Transfer Abutment - Long
Narrow Platform



NobelActive™ Compatible *	TAP-NA-NP15L	TAP-NA-NP30L	TAP-NA-NP45L
Osstem® Hiossen® Compatible *	TAP-TS-NP15L	TAP-TS-NP30L	TAP-TS-NP45L
Astra Tech™TX Compatible *	TAP-OS-SP15L	TAP-OS-SP30L	TAP-OS-SP45L

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.


Conical Connection



Transfer Abutment Standard Kit - Titanium




Angled Transfer Abutment - 15°
Regular Platform




NobelActive™ Compatible *	TAP-NA-RP1515	TAP-NA-RP1530
Osstem® Hiossen® Compatible *	TAP-TS-RP1515	TAP-TS-RP1530
Dentium® Compatible *	TAP-DE-RP1515	TAP-DE-RP1530
MegaGen® AnyOne Compatible *	TAP-AO-RP1515	TAP-AO-RP1530
MegaGen® AnyRidge Compatible *	TAP-AR-RP1515	TAP-AR-RP1530
Astra Tech™TX Compatible *	TAP-OS-WP1515	TAP-OS-WP1530
Dio Implant® Compatible *	TAP-UF-RP1515	TAP-UF-RP1530
NeoBiotech® Compatible *	TAP-IS-RP1515	TAP-IS-RP1530

Angled Transfer Abutment - 25°
Regular Platform




NobelActive™ Compatible *	TAP-NA-RP2515	TAP-NA-RP2530
Osstem® Hiossen® Compatible *	TAP-TS-RP2515	TAP-TS-RP2530
Dentium® Compatible *	TAP-DE-RP2515	TAP-DE-RP2530
MegaGen® AnyOne Compatible *	TAP-AO-RP2515	TAP-AO-RP2530
MegaGen® AnyRidge Compatible *	TAP-AR-RP2515	TAP-AR-RP2530
Astra Tech™TX Compatible *	TAP-OS-WP2515	TAP-OS-WP2530
Dio Implant® Compatible *	TAP-UF-RP2515	TAP-UF-RP2530
NeoBiotech® Compatible *	TAP-IS-RP2515	TAP-IS-RP2530

Angled Transfer Abutment - 15°
Narrow Platform




NobelActive™ Compatible *	TAP-NA-NP1515	TAP-NA-NP1530
Osstem® Hiossen® Compatible *	TAP-TS-NP1515	TAP-TS-NP1530
Astra Tech™TX Compatible *	TAP-OS-SP1515	TAP-OS-SP1530

Angled Transfer Abutment - 25°
Narrow Platform



NobelActive™ Compatible *	TAP-NA-NP2515	TAP-NA-NP2530
Osstem® Hiossen® Compatible *	TAP-TS-NP2515	TAP-TS-NP2530
Astra Tech™TX Compatible *	TAP-OS-SP2515	TAP-OS-SP2530

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.


 Professional Kit
Regular Platform

NobelActive™ Compatible *	LOCW-NA-RPPK
Astra Tech™TX Compatible *	LOCW-OS-WPPK
Dio Implant® Compatible *	LOCW-UF-RPPK
NeoBiotech® Compatible *	LOCW-IS-RPPK

Osstem® Hiossen® Compatible *	LOCW-TS-RPPK
Dentium® Compatible *	LOCW-DE-RPPK
MegaGen® AnyOne Compatible *	LOCW-AO-RPPK

 Professional Kit
Narrow Platform

NobelActive™ Compatible *	LOCW-NA-NPPK
Osstem® Hiossen® Compatible *	LOCW-TS-NPPK
Astra Tech™TX Compatible *	LOCW-OS-SPPK




Pop Clicq™ - Components

 LOCW-RCV Violet 2.7kg (4 pcs.)	 LOCW-RCT Transparent 1.8kg (4 pcs.)	 LOCW-RCP Pink 1.2kg (4 pcs.)	 LOCW-RCY Yellow 0.6kg (4 pcs.)	 LOCW-RCB Black laboratory (4 pcs.)
Metal Housing  LOCW-MH	Pick-Up Impression Transfer Assembled with Yellow Retentive Cap  LOCW-PU	Analog for Abutment  LOCW-AN	Protective Disc  LOCW-PD	
				OD-TK Overdenture Tool KIT for Insertion & Extraction of Retentive Caps Suitable for all Overdenture lines



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



 Starter Kit
Regular Platform

NobelActive™ Compatible *	LOCW-NA-RPSK
Astra Tech™TX Compatible *	LOCW-OS-WPSK
Dio Implant® Compatible *	LOCW-UF-RPSK
NeoBiotech® Compatible *	LOCW-IS-RPSK

Osstem® Hiossen® Compatible *	LOCW-TS-RPSK
Dentium® Compatible *	LOCW-DE-RPSK
MegaGen® AnyOne Compatible *	LOCW-AO-RPSK

 Starter Kit
Narrow Platform

NobelActive™ Compatible *	LOCW-NA-NPSK
Osstem® Hiossen® Compatible *	LOCW-TS-NPSK
Astra Tech™TX Compatible *	LOCW-OS-SPSK



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Pop Clicq™
Straight Mini Kit



Straight Titanium Pop-Clicq™ - Mini Kit Regular Platform

H	1mm	2mm	3mm	4mm	5mm	6mm
NobelActive™ Compatible *	LOCW-NA-RP10K	LOCW-NA-RP20K	LOCW-NA-RP30K	LOCW-NA-RP40K	LOCW-NA-RP50K	LOCW-NA-RP60K
Osstem® Hiossen® Compatible *	LOCW-TS-RP10K	LOCW-TS-RP20K	LOCW-TS-RP30K	LOCW-TS-RP40K	LOCW-TS-RP50K	LOCW-TS-RP60K
Dentium® Compatible *	LOCW-DE-RP10K	LOCW-DE-RP20K	LOCW-DE-RP30K	LOCW-DE-RP40K	LOCW-DE-RP50K	LOCW-DE-RP60K
MegaGen® AnyOne Compatible *	LOCW-AO-RP10K	LOCW-AO-RP20K	LOCW-AO-RP30K	LOCW-AO-RP40K	LOCW-AO-RP50K	LOCW-AO-RP60K
MegaGen® AnyRidge Compatible *	LOCW-AR-RP10K	LOCW-AR-RP20K	LOCW-AR-RP30K	LOCW-AR-RP40K	LOCW-AR-RP50K	LOCW-AR-RP60K
Astra Tech™ TX Compatible *	LOCW-OS-WP10K	LOCW-OS-WP20K	LOCW-OS-WP30K	LOCW-OS-WP40K	LOCW-OS-WP50K	LOCW-OS-WP60K
Dio Implant® Compatible *	LOCW-UF-RP10K	LOCW-UF-RP20K	LOCW-UF-RP30K	LOCW-UF-RP40K	LOCW-UF-RP50K	LOCW-UF-RP60K
NeoBiotech® Compatible *	LOCW-IS-RP10K	LOCW-IS-RP20K	LOCW-IS-RP30K	LOCW-IS-RP40K	LOCW-IS-RP50K	LOCW-IS-RP60K



Pop Clicq™ Straight Mini Kit

Straight Titanium Pop-Clicq™ - Mini Kit Narrow Platform

H	1mm	2mm	3mm	4mm	5mm	6mm
NobelActive™ Compatible *	LOCW-NA-NP10K	LOCW-NA-NP20K	LOCW-NA-NP30K	LOCW-NA-NP40K	LOCW-NA-NP50K	LOCW-NA-NP60K
Osstem® Hiossen® Compatible *	LOCW-TS-NP10K	LOCW-TS-NP20K	LOCW-TS-NP30K	LOCW-TS-NP40K	LOCW-TS-NP50K	LOCW-TS-NP60K
Astra Tech™ TX Compatible *	LOCW-OS-SP10K	LOCW-OS-SP20K	LOCW-OS-SP30K	LOCW-OS-SP40K	LOCW-OS-SP50K	LOCW-OS-SP60K

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Angled Titanium Pop-Clicq™ - Mini Kit - 18°
Regular Platform

	H	1mm	2mm	3mm
NobelActive™ Compatible*	LOCW-NA-RP1810K	LOCW-NA-RP1820K	LOCW-NA-RP1830K	
Osstem® Hiossen® Compatible*	LOCW-TS-RP1810K	LOCW-TS-RP1820K	LOCW-TS-RP1830K	
Dentium® Compatible*	LOCW-DE-RP1810K	LOCW-DE-RP1820K	LOCW-DE-RP1830K	
MegaGen® AnyOne Compatible*	LOCW-AO-RP1810K	LOCW-AO-RP1820K	LOCW-AO-RP1830K	
MegaGen® AnyRidge Compatible*	LOCW-AR-RP1810K	LOCW-AR-RP1820K	LOCW-AR-RP1830K	
Astra Tech™TX Compatible*	LOCW-OS-WP1810K	LOCW-OS-WP1820K	LOCW-OS-WP1830K	
Dio Implant® Compatible*	LOCW-UF-RP1810K	LOCW-UF-RP1820K	LOCW-UF-RP1830K	
NeoBiotech® Compatible*	LOCW-IS-RP1810K	LOCW-IS-RP1820K	LOCW-IS-RP1830K	



**Pop Clicq™
Angled Mini Kit**

Angled Titanium Pop-Clicq™ - Mini Kit - 18°
Narrow Platform

	H	1mm	2mm	3mm
NobelActive™ Compatible*	LOCW-NA-NP1810K	LOCW-NA-NP1820K	LOCW-NA-NP1830K	
Osstem® Hiossen® Compatible*	LOCW-TS-NP1810K	LOCW-TS-NP1820K	LOCW-TS-NP1830K	
Astra Tech™TX Compatible*	LOCW-OS-SP1810K	LOCW-OS-SP1820K	LOCW-OS-SP1830K	



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.




Angled Titanium Pop-Clicq™ - Mini Kit - 30°
Regular Platform



	1mm	2mm	3mm
H			
NobelActive™ Compatible *	LOCW-NA-RP3010K	LOCW-NA-RP3020K	LOCW-NA-RP3030K
Osstem® Hiossen® Compatible *	LOCW-TS-RP3010K	LOCW-TS-RP3020K	LOCW-TS-RP3030K
Dentium® Compatible *	LOCW-DE-RP3010K	LOCW-DE-RP3020K	LOCW-DE-RP3030K
MegaGen® AnyOne Compatible *	LOCW-AO-RP3010K	LOCW-AO-RP3020K	LOCW-AO-RP3030K
MegaGen® AnyRidge Compatible *	LOCW-AR-RP3010K	LOCW-AR-RP3020K	LOCW-AR-RP3030K
Astra Tech™TX Compatible *	LOCW-OS-WP3010K	LOCW-OS-WP3020K	LOCW-OS-WP3030K
Dio Implant® Compatible *	LOCW-UF-RP3010K	LOCW-UF-RP3020K	LOCW-UF-RP3030K
NeoBiotech® Compatible *	LOCW-IS-RP3010K	LOCW-IS-RP3020K	LOCW-IS-RP3030K



Angled Titanium Pop-Clicq™ - Mini Kit - 30°
Narrow Platform



	1mm	2mm	3mm
H			
NobelActive™ Compatible *	LOCW-NA-NP3010K	LOCW-NA-NP3020K	LOCW-NA-NP3030K
Osstem® Hiossen® Compatible *	LOCW-TS-NP3010K	LOCW-TS-NP3020K	LOCW-TS-NP3030K
Astra Tech™TX Compatible *	LOCW-OS-SP3010K	LOCW-OS-SP3020K	LOCW-OS-SP3030K

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Pop Clicq™ - Bar



Design Patent Pending No. 62154



Get a permanent bar for your patient
within minutes



Compatible with:



Conical Connection



Internal Hex

Conical Connection



Pop Clicq™ - Bar



Kits REF (Type A)

LOCW-BK2

LOCW-BK2D (Screw Driver Included)

Kits REF (Type B)

LOCW-BK4

LOCW-BK4D (Screw Driver Included)

Pop Clicq™ - Bar Components



Bars *2



(A) Adjustable Metal Housing *2
(B) Adjustable Metal Housing *4
(With a fixed screw)



Pop Clicq TM Bar
(A) Metal Housing Tool \ Adapter *2
(B) Metal Housing Tool \ Adapter *4



(A) Analog for Pop Clicq™ Abutment *2
(B) Analog for Pop Clicq™ Abutment *4



(A) Pick-Up Impression Transfer *2
(B) Pick-Up Impression Transfer *4



Adjusting Screws for
(A) Metal Housing *3 (Incl. spares)
(B) Metal Housing *6 (Incl. spares)

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.




Professional Kit
Regular Platform

NobelActive™ Compatible *	MUL-NA-RPPK
Astra Tech™TX Compatible *	MUL-OS-WPPK
Dio Implant® Compatible *	MUL-UF-RPPK
NeoBiotech® Compatible *	MUL-IS-RPPK


Osstem® Hiossen® Compatible *	MUL-TS-RPPK
Dentium® Compatible *	MUL-DE-RPPK
MegaGen® AnyOne Compatible *	MUL-AO-RPPK

Professional Kit
Narrow Platform

NobelActive™ Compatible *	MUL-NA-NPPK
Osstem® Hiossen® Compatible *	MUL-TS-NPPK
Astra Tech™TX Compatible *	MUL-OS-SPPK



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Starter Kit
Regular Platform

NobelActive™ Compatible *	MUL-NA-RPSK
Astra Tech™TX Compatible *	MUL-OS-WPSK
Dio Implant® Compatible *	MUL-UF-RPSK
NeoBiotech® Compatible *	MUL-IS-RPSK

Osstem® Hiossen® Compatible *	MUL-TS-RPSK
Dentium® Compatible *	MUL-DE-RPSK
MegaGen® AnyOne Compatible *	MUL-AO-RPSK

Starter Kit
Narrow Platform

NobelActive™ Compatible *	MUL-NA-NPSK
Osstem® Hiossen® Compatible *	MUL-TS-NPSK
Astra Tech™TX Compatible *	MUL-OS-SPSK



Optional:
* Rotational Adhesive cap
* Anti-Rotational Adhesive cap



Abutment Analog



MUL-KS-AN

Closed Tray Plastic Transfer



MUL-KS-PUP
[MUL-KS-TSSV - included]
MUL-KS-PUPUG
[MUL-KS-TSUG - included]
MUL-KS-PUPAR
[MUL-KS-TSAR - included]

Castable Plastic Sleeve



MUL-KS-SLVP
[MUL-KS-TSSV - included]
MUL-KS-SLVPUG
[MUL-KS-TSUG - included]
MUL-KS-SLVPAR
[MUL-KS-TSAR - included]

Cobalt Sleeve UCLA



MUL-KS-SLVCC
[MUL-KS-TSSV - included]
MUL-KS-SLVCCUG
[MUL-KS-TSUG - included]
MUL-KS-SLVCCR
[MUL-KS-TSAR - included]

Open Tray Transfer



MUL-KS-OT
[MUL-KS-SOT - included]
MUL-KS-OTUG
[MUL-KS-SOTUG - included]
MUL-KS-OTAR
[MUL-KS-SOTAR - included]

Healing cap



MUL-KS-HA
[MUL-KS-TSSV - included]
MUL-KS-HAUG
[MUL-KS-TSUG - included]
MUL-KS-HAAR
[MUL-KS-TSAR - included]

Scan Body (Peek)



MUL-KS-SB
[MUL-KS-TSSV - included]
MUL-KS-SBUG
[MUL-KS-TSUG - included]
MUL-KS-SBAR
[MUL-KS-TSAR - included]

Temporary Titanium Sleeve



MUL-KS-SLVTI
[MUL-KS-TSSV - included]
MUL-KS-SLVTIUG
[MUL-KS-TSUG - included]
MUL-KS-SLVTIAR
[MUL-KS-TSAR - included]

Rotational Titanium Base



MUL-KS-TBR
[MUL-KS-TSSV - included]
MUL-KS-TBRUG
[MUL-KS-TSUG - included]
MUL-KS-TBRAR
[MUL-KS-TSAR - included]

Anti Rotational Titanium Base



MUL-KS-TB
[MUL-KS-TSSV - included]
MUL-KS-TBUG
[MUL-KS-TSUG - included]
MUL-KS-TBAR
[MUL-KS-TSAR - included]

*For Single Unit Restoration.

Rotational Adhesive Cap



MUL-KS-ADR
[MUL-KS-TSSV - included]
MUL-KS-ADRUG
[MUL-KS-TSUG - included]
MUL-KS-ADRAR
[MUL-KS-TSAR - included]

Anti Rotational Adhesive Cap



MUL-KS-AD
[MUL-KS-TSSV - included]
MUL-KS-ADUG
[MUL-KS-TSUG - included]
MUL-KS-ADAR
[MUL-KS-TSAR - included]

*For Single Unit Restoration.

Prosthetic Titanium Screw



MUL-KS-TSSV (1.28mm)
MUL-KS-TSUG (UniGrip)
MUL-KS-TSAR (1.2mm)



Conical Connection



Multi Clicq™ - On4/On6 Straight Mini Kit



Straight Mini Kit Regular Platform

H	1.5mm	2.5mm	3.5mm	4.5mm	5.5mm
NobelActive™ Compatible *	MUL-NA-RPKS15K	MUL-NA-RPKS25K	MUL-NA-RPKS35K	MUL-NA-RPKS45K	MUL-NA-RPKS55K
Osstem® Hiossen® Compatible *	MUL-TS-RPKS15K	MUL-TS-RPKS25K	MUL-TS-RPKS35K	MUL-TS-RPKS45K	MUL-TS-RPKS55K
Dentium® Compatible *	MUL-DE-RPKS15K	MUL-DE-RPKS25K	MUL-DE-RPKS35K	MUL-DE-RPKS45K	MUL-DE-RPKS55K
MegaGen® AnyOne Compatible *	MUL-AO-RPKS15K	MUL-AO-RPKS25K	MUL-AO-RPKS35K	MUL-AO-RPKS45K	MUL-AO-RPKS55K
MegaGen® AnyRidge Compatible *	MUL-AR-RPKS15K	MUL-AR-RPKS25K	MUL-AR-RPKS35K	MUL-AR-RPKS45K	MUL-AR-RPKS55K
Astra Tech™TX Compatible *	MUL-OS-WPKS15K	MUL-OS-WPKS25K	MUL-OS-WPKS35K	MUL-OS-WPKS45K	MUL-OS-WPKS55K
Dio Implant® Compatible *	MUL-UF-RPKS15K	MUL-UF-RPKS25K	MUL-UF-RPKS35K	MUL-UF-RPKS45K	MUL-UF-RPKS55K
NeoBiotech® Compatible *	MUL-IS-RPKS15K	MUL-IS-RPKS25K	MUL-IS-RPKS35K	MUL-IS-RPKS45K	MUL-IS-RPKS55K



Straight Mini Kit Narrow Platform

H	1.5mm	2.5mm	3.5mm	4.5mm	5.5mm
NobelActive™ Compatible *	MUL-NA-NPKS15K	MUL-NA-NPKS25K	MUL-NA-NPKS35K	MUL-NA-NPKS45K	MUL-NA-NPKS55K
Osstem® Hiossen® Compatible *	MUL-TS-NPKS15K	MUL-TS-NPKS25K	MUL-TS-NPKS35K	MUL-TS-NPKS45K	MUL-TS-NPKS55K
Astra Tech™TX Compatible *	MUL-OS-SPKS15K	MUL-OS-SPKS25K	MUL-OS-SPKS35K	MUL-OS-SPKS45K	MUL-OS-SPKS55K

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Multi Clicq™ - On4/On6 Angled Mini Kit



One Piece Angled 17°
Abutment Mini Kit
Regular Platform



H G	1mm 2.5mm	2mm 3.5mm
NobelActive™ Compatible *	MUL-NA-RPKS1710K	MUL-NA-RPKS1720K
Osstem® Hiossen® Compatible *	MUL-TS-RPKS1710K	MUL-TS-RPKS1720K
Dentium® Compatible *	MUL-DE-RPKS1710K	MUL-DE-RPKS1720
MegaGen® AnyOne Compatible *	MUL-AO-RPKS1710K	MUL-AO-RPKS1720K
MegaGen® AnyRidge Compatible *	MUL-AR-RPKS1710K	MUL-AR-RPKS1720K
Astra Tech™TX Compatible *	MUL-OS-WPKS1710K	MUL-OS-WPKS1720K
Dio Implant® Compatible *	MUL-UF-RPKS1710K	MUL-UF-RPKS1720K
NeoBiotech® Compatible *	MUL-IS-RPKS1710K	MUL-IS-RPKS1720K

One Piece Angled 30°
Abutment Mini Kit
Regular Platform



H G	1mm 3.5mm	2mm 4.5mm
NobelActive™ Compatible *	MUL-NA-RPKS3010K	MUL-NA-RPKS3020K
Osstem® Hiossen® Compatible *	MUL-TS-RPKS3010K	MUL-TS-RPKS3020K
Dentium® Compatible *	MUL-DE-RPKS3010K	MUL-DE-RPKS3020K
MegaGen® AnyOne Compatible *	MUL-AO-RPKS3010K	MUL-AO-RPKS3020K
MegaGen® AnyRidge Compatible *	MUL-AR-RPKS3010K	MUL-AR-RPKS3020K
Astra Tech™TX Compatible *	MUL-OS-WPKS3010K	MUL-OS-WPKS3020K
Dio Implant® Compatible *	MUL-UF-RPKS3010K	MUL-UF-RPKS3020K
NeoBiotech® Compatible *	MUL-IS-RPKS3010K	MUL-IS-RPKS3020K



One Piece Angled 17°
Abutment Mini Kit
Narrow Platform



H G	1mm 2.5mm	2mm 3.5mm
NobelActive™ Compatible *	MUL-NA-NPKS1710K	MUL-NA-NPKS1720K
Osstem® Hiossen® Compatible *	MUL-TS-NPKS1710K	MUL-TS-NPKS1720K
Astra Tech™TX Compatible *	MUL-OS-SPKS1710K	MUL-OS-SPKS1720K

One Piece Angled 30°
Abutment Mini Kit
Narrow Platform



H G	1mm 3.5mm	2mm 4.5mm
NobelActive™ Compatible *	MUL-NA-NPKS3010K	MUL-NA-NPKS3020K
Osstem® Hiossen® Compatible *	MUL-TS-NPKS3010K	MUL-TS-NPKS3020K
Astra Tech™TX Compatible *	MUL-OS-SPKS3010K	MUL-OS-SPKS3020K

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Straight Abutment
Regular Platform

	H	1.5mm	2.5mm	3.5mm	4.5mm	5.5mm
NobelActive™ Compatible*	MUL-NA-RPKS15	MUL-NA-RPKS25	MUL-NA-RPKS35	MUL-NA-RPKS45	MUL-NA-RPKS55	
Osstem® Hiossen® Compatible*	MUL-TS-RPKS15	MUL-TS-RPKS25	MUL-TS-RPKS35	MUL-TS-RPKS45	MUL-TS-RPKS55	
Dentium® Compatible*	MUL-DE-RPKS15	MUL-DE-RPKS25	MUL-DE-RPKS35	MUL-DE-RPKS45	MUL-DE-RPKS55	
MegaGen® AnyOne Compatible*	MUL-AO-RPKS15	MUL-AO-RPKS25	MUL-AO-RPKS35	MUL-AO-RPKS45	MUL-AO-RPKS55	
MegaGen® AnyRidge Compatible*	MUL-AR-RPKS15	MUL-AR-RPKS25	MUL-AR-RPKS35	MUL-AR-RPKS45	MUL-AR-RPKS55	
Astra Tech™TX Compatible*	MUL-OS-WPKS15	MUL-OS-WPKS25	MUL-OS-WPKS35	MUL-OS-WPKS45	MUL-OS-WPKS55	
Dio Implant® Compatible*	MUL-UF-RPKS15	MUL-UF-RPKS25	MUL-UF-RPKS35	MUL-UF-RPKS45	MUL-UF-RPKS55	
NeoBiotech® Compatible*	MUL-IS-RPKS15	MUL-IS-RPKS25	MUL-IS-RPKS35	MUL-IS-RPKS45	MUL-IS-RPKS55	

Straight Abutment
Narrow Platform

	H	1.5mm	2.5mm	3.5mm	4.5mm	5.5mm
NobelActive™ Compatible*	MUL-NA-NPKS15	MUL-NA-NPKS25	MUL-NA-NPKS35	MUL-NA-NPKS45	MUL-NA-NPKS55	
Osstem® Hiossen® Compatible*	MUL-TS-NPKS15	MUL-TS-NPKS25	MUL-TS-NPKS35	MUL-TS-NPKS45	MUL-TS-NPKS55	
Astra Tech™TX Compatible*	MUL-OS-SPKS15	MUL-OS-SPKS25	MUL-OS-SPKS35	MUL-OS-SPKS45	MUL-OS-SPKS55	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Multi Clicq™ - On4/On6 Standard Kit



One Piece Angled Abutment - 17° Regular Platform

H		1mm 2.5mm	2mm 3.5mm
G			
NobelActive™ Compatible *	MUL-NA-RPKS1710	MUL-NA-RPKS1720	
Osstem® Hiossen® Compatible *	MUL-TS-RPKS1710	MUL-TS-RPKS1720	
Dentium® Compatible *	MUL-DE-RPKS1710	MUL-DE-RPKS1720	
MegaGen® AnyOne Compatible *	MUL-AO-RPKS1710	MUL-AO-RPKS1720	
MegaGen® AnyRidge Compatible *	MUL-AR-RPKS1710	MUL-AR-RPKS1720	
Astra Tech™TX Compatible *	MUL-OS-WPKS1710	MUL-OS-WPKS1720	
Dio Implant® Compatible *	MUL-UF-RPKS1710	MUL-UF-RPKS1720	
NeoBiotech® Compatible *	MUL-IS-RPKS1710	MUL-IS-RPKS1720	

One Piece Angled Abutment - 30° Regular Platform

H		1mm 3.5mm	2mm 4.5mm
G			
NobelActive™ Compatible *	MUL-NA-RPKS3010	MUL-NA-RPKS3020	
Osstem® Hiossen® Compatible *	MUL-TS-RPKS3010	MUL-TS-RPKS3020	
Dentium® Compatible *	MUL-DE-RPKS3010	MUL-DE-RPKS3020	
MegaGen® AnyOne Compatible *	MUL-AO-RPKS3010	MUL-AO-RPKS3020	
MegaGen® AnyRidge Compatible *	MUL-AR-RPKS3010	MUL-AR-RPKS3020	
Astra Tech™TX Compatible *	MUL-OS-WPKS3010	MUL-OS-WPKS3020	
Dio Implant® Compatible *	MUL-UF-RPKS3010	MUL-UF-RPKS3020	
NeoBiotech® Compatible *	MUL-IS-RPKS3010	MUL-IS-RPKS3020	

One Piece Angled Abutment - 17° Narrow Platform

H		1mm 2.5mm	2mm 3.5mm
G			
NobelActive™ Compatible *	MUL-NA-NPKS1710	MUL-NA-NPKS1720	
Osstem® Hiossen® Compatible *	MUL-TS-NPKS1710	MUL-TS-NPKS1720	
Astra Tech™TX Compatible *	MUL-OS-SPKS1710	MUL-OS-SPKS1720	

One Piece Angled Abutment - 30° Narrow Platform

H		1mm 3.5mm	2mm 4.5mm
G			
NobelActive™ Compatible *	MUL-NA-NPKS3010	MUL-NA-NPKS3020	
Osstem® Hiossen® Compatible *	MUL-TS-NPKS3010	MUL-TS-NPKS3020	
Astra Tech™TX Compatible *	MUL-OS-SPKS3010	MUL-OS-SPKS3020	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Multi Clicq™ - On4/On6 Standard Kit - Anti-Rotational



One Piece Angled Abutment - 17° Regular Platform

	H	G	H	G	1mm 2.5mm	2mm 3.5mm
NobelActive™ Compatible *			MUL-NA-RPKS1710H	MUL-NA-RPKS1720H		
Osstem® Hiossen® Compatible *			MUL-TS-RPKS1710H	MUL-TS-RPKS1720H		
Dentium® Compatible *			MUL-DE-RPKS1710H	MUL-DE-RPKS1720H		
MegaGen® AnyOne Compatible *			MUL-AO-RPKS1710H	MUL-AO-RPKS1720H		
MegaGen® AnyRidge Compatible *			MUL-AR-RPKS1710H	MUL-AR-RPKS1720H		
Astra Tech™TX Compatible *			MUL-OS-WPKS1710H	MUL-OS-WPKS1720H		
Dio Implant® Compatible *			MUL-UF-RPKS1710H	MUL-UF-RPKS1720H		
NeoBiotech® Compatible *			MUL-IS-RPKS1710H	MUL-IS-RPKS1720H		

One Piece Angled Abutment - 30° Regular Platform

	H	G	H	G	1mm 3.5mm	2mm 4.5mm
NobelActive™ Compatible *			MUL-NA-RPKS3010H	MUL-NA-RPKS3020H		
Osstem® Hiossen® Compatible *			MUL-TS-RPKS3010H	MUL-TS-RPKS3020H		
Dentium® Compatible *			MUL-DE-RPKS3010H	MUL-DE-RPKS3020H		
MegaGen® AnyOne Compatible *			MUL-AO-RPKS3010H	MUL-AO-RPKS3020H		
MegaGen® AnyRidge Compatible *			MUL-AR-RPKS3010H	MUL-AR-RPKS3020H		
Astra Tech™TX Compatible *			MUL-OS-WPKS3010H	MUL-OS-WPKS3020H		
Dio Implant® Compatible *			MUL-UF-RPKS3010H	MUL-UF-RPKS3020H		
NeoBiotech® Compatible *			MUL-IS-RPKS3010H	MUL-IS-RPKS3020H		

One Piece Angled Abutment - 17° Narrow Platform

	H	G	H	G	1mm 2.5mm	2mm 3.5mm
NobelActive™ Compatible *			MUL-NA-NPKS1710H	MUL-NA-NPKS1720H		
Osstem® Hiossen® Compatible *			MUL-TS-NPKS1710H	MUL-TS-NPKS1720H		
Astra Tech™TX Compatible *			MUL-OS-SPKS1710H	MUL-OS-SPKS1720H		

One Piece Angled Abutment - 30° Narrow Platform


	H	G	H	G	1mm 3.5mm	2mm 4.5mm
NobelActive™ Compatible *			MUL-NA-NPKS3010H	MUL-NA-NPKS3020H		
Osstem® Hiossen® Compatible *			MUL-TS-NPKS3010H	MUL-TS-NPKS3020H		
Astra Tech™TX Compatible *			MUL-OS-SPKS3010H	MUL-OS-SPKS3020H		

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Clicq™ - Multi Bar



Get a permanent bar
for your patient
within minutes



Conical Connection

Design Patent Pending No. 62154
Patent Pending No. US 62/679,389



Clicq™ - Multi Bar
On4/On6 Kit



Clicq™ - Multi Bar On4 Kit
Regular Platform



(Screw Driver Included)

NobelActive™ Compatible *	MUL-NA-BK4	MUL-NA-BK4D
Osstem® Hiossen® Compatible *	MUL-TS-BK4	MUL-TS-BK4D
Dentium® Compatible *	MUL-DE-BK4	MUL-DE-BK4D
MegaGen® AnyOne Compatible *	MUL-AO-BK4	MUL-AO-BK4D
MegaGen® AnyRidge Compatible *	MUL-AR-BK4	MUL-AR-BK4D
Astra Tech™TX Compatible *	MUL-OS-BK4	MUL-OS-BK4D
Dio Implant® Compatible *	MUL-UF-BK4	MUL-UF-BK4D
NeoBiotech® Compatible *	MUL-IS-BK4	MUL-IS-BK4D

Clicq™ - Multi Bar On6 Kit
Regular Platform



(Screw Driver Included)

NobelActive™ Compatible *	MUL-NA-BK6	MUL-NA-BK6D
Osstem® Hiossen® Compatible *	MUL-TS-BK6	MUL-TS-BK6D
Dentium® Compatible *	MUL-DE-BK6	MUL-DE-BK6D
MegaGen® AnyOne Compatible *	MUL-AO-BK6	MUL-AO-BK6D
MegaGen® AnyRidge Compatible *	MUL-AR-BK6	MUL-AR-BK6D
Astra Tech™TX Compatible *	MUL-OS-BK6	MUL-OS-BK6D
Dio Implant® Compatible *	MUL-UF-BK6	MUL-UF-BK6D
NeoBiotech® Compatible *	MUL-IS-BK6	MUL-IS-BK6D

Clicq™ - Multi Bar Components



Bars *2



Adjusting Micro Screw *2
(spares)



Prosthetic KS Ti-Screw
On 4 - *6 (Incl. 2 spares)
On 6 - *8 (Incl. 2 spares)



Pre-Fabricated Permanent
Titanium Adjustable Sleeve *4/6
W/ Adjusting Micro Screw *8/12




Hand Screw Driver

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.




 Straight Anatomic Abutment
Regular Platform



	1.5mm 9.5mm	3mm 11mm	4.5mm 12.5mm
NobelActive™ Compatible *	ZAA-NA-RP15	ZAA-NA-RP30	ZAA-NA-RP45
Osstem® Hiossen® Compatible *	ZAA-TS-RP15	ZAA-TS-RP30	ZAA-TS-RP45
Dentium® Compatible *	ZAA-DE-RP15	ZAA-DE-RP30	ZAA-DE-RP45
MegaGen® AnyOne Compatible *	ZAA-AO-RP15	ZAA-AO-RP30	ZAA-AO-RP45
MegaGen® AnyRidge Compatible *	ZAA-AR-RP15	ZAA-AR-RP30	ZAA-AR-RP45
Astra Tech™TX Compatible *	ZAA-OS-WP15	ZAA-OS-WP30	ZAA-OS-WP45
Dio Implant® Compatible *	ZAA-UF-RP15	ZAA-UF-RP30	ZAA-UF-RP45
NeoBiotech® Compatible *	ZAA-IS-RP15	ZAA-IS-RP30	ZAA-IS-RP45

 Straight Anatomic Abutment
Narrow Platform



	1.5mm 9.5mm	3mm 11mm	4.5mm 12.5mm
NobelActive™ Compatible *	ZAA-NA-NP15	ZAA-NA-NP30	ZAA-NA-NP45
Osstem® Hiossen® Compatible *	ZAA-TS-NP15	ZAA-TS-NP30	ZAA-TS-NP45
Astra Tech™TX Compatible *	ZAA-OS-SP15	ZAA-OS-SP30	ZAA-OS-SP45

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Angled Anatomic Abutment - 15°
Regular Platform

H	1.5mm 8.5mm	3mm 10mm	4.5mm 11.5mm
NobelActive™ Compatible *	ZAA-NA-RP1515 ZAA-NA-RP1530 ZAA-NA-RP1545		
Osstem® Hiossen® Compatible *	ZAA-TS-RP1515 ZAA-TS-RP1530 ZAA-TS-RP1545		
Dentium® Compatible *	ZAA-DE-RP1515 ZAA-DE-RP1530 ZAA-DE-RP1545		
MegaGen® AnyOne Compatible *	ZAA-AO-RP1515 ZAA-AO-RP1530 ZAA-AO-RP1545		
MegaGen® AnyRidge Compatible *	ZAA-AR-RP1515 ZAA-AR-RP1530 ZAA-AR-RP1545		
Astra Tech™TX Compatible *	ZAA-OS-WP1515 ZAA-OS-WP1530 ZAA-OS-WP1545		
Dio Implant® Compatible *	ZAA-UF-RP1515 ZAA-UF-RP1530 ZAA-UF-RP1545		
NeoBiotech® Compatible *	ZAA-IS-RP1515 ZAA-IS-RP1530 ZAA-IS-RP1545		

Angled Anatomic Abutment - 25°
Regular Platform

H	1.5mm 8.5mm	3mm 10mm	4.5mm 11.5mm
NobelActive™ Compatible *	ZAA-NA-RP2515 ZAA-NA-RP2530 ZAA-NA-RP2545		
Osstem® Hiossen® Compatible *	ZAA-TS-RP2515 ZAA-TS-RP2530 ZAA-TS-RP2545		
Dentium® Compatible *	ZAA-DE-RP2515 ZAA-DE-RP2530 ZAA-DE-RP2545		
MegaGen® AnyOne Compatible *	ZAA-AO-RP2515 ZAA-AO-RP2530 ZAA-AO-RP2545		
MegaGen® AnyRidge Compatible *	ZAA-AR-RP2515 ZAA-AR-RP2530 ZAA-AR-RP2545		
Astra Tech™TX Compatible *	ZAA-OS-WP2515 ZAA-OS-WP2530 ZAA-OS-WP2545		
Dio Implant® Compatible *	ZAA-UF-RP2515 ZAA-UF-RP2530 ZAA-UF-RP2545		
NeoBiotech® Compatible *	ZAA-IS-RP2515 ZAA-IS-RP2530 ZAA-IS-RP2545		

Angled Anatomic Abutment - 15°
Narrow Platform

H	1.5mm 8.5mm	3mm 10mm	4.5mm 11.5mm
NobelActive™ Compatible *	ZAA-NA-NP1515 ZAA-NA-NP1530 ZAA-NA-NP1545		
Osstem® Hiossen® Compatible *	ZAA-TS-NP1515 ZAA-TS-NP1530 ZAA-TS-NP1545		
Astra Tech™TX Compatible *	ZAA-OS-SP1515 ZAA-OS-SP1530 ZAA-OS-SP1545		

Angled Anatomic Abutment - 25°
Narrow Platform

H	1.5mm 8.5mm	3mm 10mm	4.5mm 11.5mm
NobelActive™ Compatible *	ZAA-NA-NP2515 ZAA-NA-NP2530 ZAA-NA-NP2545		
Osstem® Hiossen® Compatible *	ZAA-TS-NP2515 ZAA-TS-NP2530 ZAA-TS-NP2545		
Astra Tech™TX Compatible *	ZAA-OS-SP2515 ZAA-OS-SP2530 ZAA-OS-SP2545		

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Straight Anatomic Abutment
Regular Platform

H		
L	1.5mm	9.5mm
NobelActive™ Compatible*	ZAAU-NA-RP15	
Osstem® Hiossen® Compatible*	ZAAU-TS-RP15	
Dentium® Compatible*	ZAAU-DE-RP15	
MegaGen® AnyOne Compatible*	ZAAU-AO-RP15	
MegaGen® AnyRidge Compatible*	ZAAU-AR-RP15	
Astra Tech™TX Compatible*	ZAAU-OS-WP15	
Dio Implant® Compatible*	ZAAU-UF-RP15	
NeoBiotech® Compatible*	ZAAU-IS-RP15	

Angled Anatomic Abutment - 15°
Regular Platform

H		
L	1.5mm	9.5mm
NobelActive™ Compatible*	ZAAU-NA-RP1515	
Osstem® Hiossen® Compatible*	ZAAU-TS-RP1515	
Dentium® Compatible*	ZAAU-DE-RP1515	
MegaGen® AnyOne Compatible*	ZAAU-AO-RP1515	
MegaGen® AnyRidge Compatible*	ZAAU-AR-RP1515	
Astra Tech™TX Compatible*	ZAAU-OS-WP1515	
Dio Implant® Compatible*	ZAAU-UF-RP1515	
NeoBiotech® Compatible*	ZAAU-IS-RP1515	

Straight Anatomic Abutment
Narrow Platform

H		
L	1.5mm	9.5mm
NobelActive™ Compatible*	ZAAU-NA-NP15	
Osstem® Hiossen® Compatible*	ZAAU-TS-NP15	
Astra Tech™TX Compatible*	ZAAU-OS-SP15	

Angled Anatomic Abutment - 15°
Narrow Platform

H		
L	1.5mm	9.5mm
NobelActive™ Compatible*	ZAAU-NA-NP1515	
Osstem® Hiossen® Compatible*	ZAAU-TS-NP1515	
Astra Tech™TX Compatible*	ZAAU-OS-SP1515	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Titanium Abutment
Regular Platform

H	10.5mm	10.5mm (Castable Plastic Sleeve included)	12mm	12mm	L 11
NobelActive™ Compatible *	ABZ-NA-RPTI	ABZ-NA-RPTIR			
Osstem® Hiossen® Compatible *	ABZ-TS-RPTI	ABZ-TS-RPTIR			
Dentium® Compatible *	ABZ-DE-RPTI	ABZ-DE-RPTIR			
MegaGen® AnyOne Compatible *	ABZ-AO-RPTI	ABZ-AO-RPTIR			
MegaGen® AnyRidge Compatible *			ABZ-AR-RPTI	ABZ-AR-RPTIR	
Astra Tech™TX Compatible *	ABZ-OS-WPTI	ABZ-OS-WPTIR			
Dio Implant® Compatible *	ABZ-UF-RPTI	ABZ-UF-RPTIR			
NeoBiotech® Compatible *	ABZ-IS-RPTI	ABZ-IS-RPTIR			

(Castable Plastic Sleeve included)

Titanium Abutment
Narrow Platform

H	10.5mm	10.5mm (Castable Plastic Sleeve included)
NobelActive™ Compatible *	ABZ-NA-NPTI	ABZ-NA-NPTIR
Osstem® Hiossen® Compatible *	ABZ-TS-NPTI	ABZ-TS-NPTIR
Astra Tech™TX Compatible *	ABZ-OS-SPTI	ABZ-OS-SPTIR

(Castable Plastic Sleeve included)

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Chrome Cobalt Abutment
Regular Platform

H	10.5mm	10.5mm (Castable Plastic Sleeve included)	12mm	12mm	L 11
NobelActive™ Compatible *	ABZ-NA-RPCC	ABZ-NA-RPCCR			
Osstem® Hiossen® Compatible *	ABZ-TS-RPCC	ABZ-TS-RPCCR			
Dentium® Compatible *	ABZ-DE-RPCC	ABZ-DE-RPCCR			
MegaGen® AnyOne Compatible *	ABZ-AO-RPCC	ABZ-AO-RPCCR			
MegaGen® AnyRidge Compatible *			ABZ-AR-RPCC	ABZ-AR-RPCCR	
Astra Tech™TX Compatible *	ABZ-OS-WPCC	ABZ-OS-WPCCR			
Dio Implant® Compatible *	ABZ-UF-RPCC	ABZ-UF-RPCCR			
NeoBiotech® Compatible *	ABZ-IS-RPCC	ABZ-IS-RPCCR			

(Castable Plastic Sleeve included)

Chrome Cobalt Abutment
Narrow Platform

H	10.5mm	10.5mm (Castable Plastic Sleeve included)
NobelActive™ Compatible *	ABZ-NA-NPCC	ABZ-NA-NPCCR
Osstem® Hiossen® Compatible *	ABZ-TS-NPCC	ABZ-TS-NPCCR
Astra Tech™TX Compatible *	ABZ-OS-SPCC	ABZ-OS-SPCCR

(Castable Plastic Sleeve included)


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection




Casting Abutment Plastic

Plastic Abutment for Casting Regular Platform



NobelActive™ Compatible*	PAC-NA-RP	PAC-NA-RPR
Osstem® Hiossen® Compatible*	PAC-TS-RP	PAC-TS-RPR
Dentium® Compatible*	PAC-DE-RP	PAC-DE-RPR
MegaGen® AnyOne Compatible *	PAC-AO-RP	PAC-AO-RPR
MegaGen® AnyRidge Compatible *	PAC-AR-RP	PAC-AR-RPR
Astra Tech™TX Compatible *	PAC-OS-WP	PAC-OS-WPR
Dio Implant® Compatible *	PAC-UF-RP	PAC-UF-RPR
NeoBiotech® Compatible *	PAC-IS-RP	PAC-IS-RPR

Plastic Abutment for Casting Narrow Platform



NobelActive™ Compatible*	PAC-NA-NP	PAC-NA-NPR
Osstem® Hiossen® Compatible*	PAC-TS-NP	PAC-TS-NPR
Astra Tech™TX Compatible *	PAC-OS-SP	PAC-OS-SPR

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.




PEEK Scan Abutment
Regular Platform

NobelActive™ Compatible *		SB-NA-RP			
Osstem® Hiossen® Compatible *			SB-TS-RP		
Dentium® Compatible *				SB-DE-RP	
MegaGen® AnyOne Compatible *					SB-AO-RP
MegaGen® AnyRidge Compatible *					SB-AR-RP
Astra Tech™TX Compatible *		SB-OS-WP			
Dio Implant® Compatible *		SB-UF-RP			
NeoBiotech® Compatible *		SB-IS-RP			

PEEK Scan Abutment
Narrow Platform

NobelActive™ Compatible *		SB-NA-NP
Osstem® Hiossen® Compatible *		SB-TS-NP
Astra Tech™TX Compatible *		SB-OS-SP

Scan Body (Non PEEK)
for Regular & Narrow Platform



Sirona CEREC®
Compatible *



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Titanium Base
Regular Platform

		Sirona CEREC® Compatible*				
H	Ø	1.5mm Ø 4.5mm	0.8mm Ø 4.5mm	0.6mm Ø 4mm	1.5mm Ø 4.5mm	
NobelActive™ Compatible*		TB-NA-RP	TBC-NA-RP			
Osstem® Hiossen® Compatible*			TBC-TS-RP			
Dentium® Compatible*			TBC-DE-RP			
MegaGen® AnyOne Compatible*			TBC-AO-RP	TB-AO-RP4045	TB-AO-RP4545	
MegaGen® AnyRidge Compatible*			TBC-AR-RP	TB-AR-RP4045	TB-AR-RP4545	
Astra Tech™TX Compatible*			TBC-OS-WP			
Dio Implant® Compatible*			TBC-UF-RP			
NeoBiotech® Compatible*			TBC-IS-RP			

Angled Titanium Base 10°
Regular Platform

		Sirona CEREC® Compatible*				
L		10°	H			
				0.5mm	0.5mm	
NobelActive™ Compatible*		TBC-NA-RP1000	TBC-NA-RP1000R			
Osstem® Hiossen® Compatible*		TBC-TS-RP1000	TBC-TS-RP1000R			
Dentium® Compatible*		TBC-DE-RP1000	TBC-DE-RP1000R			
MegaGen® AnyOne Compatible*		TBC-AO-RP1000	TBC-AO-RP1000R			
MegaGen® AnyRidge Compatible*		TBC-AR-RP1000	TBC-AR-RP1000R			
Astra Tech™TX Compatible*		TBC-OS-WP1000	TBC-OS-WP1000R			
Dio Implant® Compatible*		TBC-UF-RP1000	TBC-UF-RP1000R			
NeoBiotech® Compatible*		TBC-IS-RP1000	TBC-IS-RP1000R			

Titanium Base
Narrow Platform

		Sirona CEREC® Compatible*				
H	Ø	1.5mm Ø 4.5mm	0.8mm Ø 4.5mm			
NobelActive™ Compatible*		TB-NA-NP	TBC-NA-NP			
Osstem® Hiossen® Compatible*			TBC-TS-NP			
Astra Tech™TX Compatible*			TBC-OS-SP			

Angled Titanium Base 10°
Narrow Platform

		Sirona CEREC® Compatible*				
L		10°	H			
				0.6mm	0.6mm	
NobelActive™ Compatible*		TBC-NA-NP1000	TBC-NA-NP1000R			
Osstem® Hiossen® Compatible*		TBC-TS-NP1000	TBC-NA-NP1000R			
Astra Tech™TX Compatible*		TBC-OS-SP1000	TBC-NA-SP1000R			

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Pre-Milled Titanium Blank
Regular Platform

	Arum Compatible *	Medentika Compatible *	NT-Trading Compatible *
NobelActive™ Compatible *	PMAM-NA-RP	PMMD-NA-RP	PMNT-NA-RP
Osstem® Hiossen® Compatible *	PMAM-TS-RP	PMMD-TS-RP	PMNT-TS-RP
Dentium® Compatible *	PMAM-DE-RP	PMMD-DE-RP	PMNT-DE-RP
MegaGen® AnyOne Compatible *	PMAM-AO-RP	PMMD-AO-RP	PMNT-AO-RP
MegaGen® AnyRidge Compatible *	PMAM-AR-RP	PMMD-AR-RP	PMNT-AR-RP
Astra Tech™TX Compatible *	PMAM-OS-WP	PMMD-OS-WP	PMNT-OS-WP
Dio Implant® Compatible *	PMAM-UF-RP	PMMD-UF-RP	PMNT-UF-RP
NeoBiotech® Compatible *	PMAM-IS-RP	PMMD-IS-RP	PMNT-IS-RP

Pre-Milled Titanium Blank
Narrow Platform

	Arum Compatible *	Medentika Compatible *	NT-Trading Compatible *
NobelActive™ Compatible *	PMAM-NA-NP	PMMD-NA-NP	PMNT-NA-NP
Osstem® Hiossen® Compatible *	PMAM-TS-NP	PMMD-TS-NP	PMNT-TS-NP
Astra Tech™TX Compatible *	PMAM-OS-SP	PMMD-OS-SP	PMNT-OS-SP

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



 Pre-Milled Titanium Blank - Rotational
Regular Platform

	Arum Compatible *	Medentika Compatible *	NT-Trading Compatible *
NobelActive™ Compatible *	PMAM-NA-RPR	PMMD-NA-RPR	PMNT-NA-RPR
Osstem® Hiossen® Compatible *	PMAM-TS-RPR	PMMD-TS-RPR	PMNT-TS-RPR
Dentium® Compatible *	PMAM-DE-RPR	PMMD-DE-RPR	PMNT-DE-RPR
MegaGen® AnyOne Compatible *	PMAM-AO-RPR	PMMD-AO-RPR	PMNT-AO-RPR
MegaGen® AnyRidge Compatible *	PMAM-AR-RPR	PMMD-AR-RPR	PMNT-AR-RPR
Astra Tech™TX Compatible *	PMAM-OS-WPR	PMMD-OS-WPR	PMNT-OS-WPR
Dio Implant® Compatible *	PMAM-UF-RPR	PMMD-UF-RPR	PMNT-UF-RPR
NeoBiotech® Compatible *	PMAM-IS-RPR	PMMD-IS-RPR	PMNT-IS-RPR

 Pre-Milled Titanium Blank - Rotational
Narrow Platform


	Arum Compatible *	Medentika Compatible *	NT-Trading Compatible *
NobelActive™ Compatible *	PMAM-NA-NPR	PMMD-NA-NPR	PMNT-NA-NPR
Osstem® Hiossen® Compatible *	PMAM-TS-NPR	PMMD-TS-NPR	PMNT-TS-NPR
Astra Tech™TX Compatible *	PMAM-OS-SPR	PMMD-OS-SPR	PMNT-OS-SPR

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection




Abutment Gauge Mini Kit



Abutment Gauge Kit Regular Platform

NobelActive™ Compatible *	AG-NA-RPK
Osstem® Hiossen® Compatible *	AG-TS-RPK
Dentium® Compatible *	AG-DE-RPK
MegaGen® AnyOne Compatible *	AG-AO-RPK
MegaGen® AnyRidge Compatible *	AG-AR-RPK
Astra Tech™TX Compatible *	AG-OS-WPK
Dio Implant® Compatible *	AG-UF-RPK
NeoBiotech® Compatible *	AG-IS-RPK



Abutment Gauge Kit Narrow Platform

NobelActive™ Compatible *	AG-NA-NPK
Osstem® Hiossen® Compatible *	AG-TS-NPK
Astra Tech™TX Compatible *	AG-OS-SPK

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Conical Connection



Abutment Gauge

Abutment Gauge Straight Regular Platform

	H	1-5mm
NobelActive™ Compatible *	AG-NA-RP00	
Osstem® Hiossen® Compatible *	AG-TS-RP00	
Dentium® Compatible *	AG-DE-RP00	
MegaGen® AnyOne Compatible *	AG-AO-RP00	
MegaGen® AnyRidge Compatible *	AG-AR-RP00	
Astra Tech™TX Compatible *	AG-OS-WP00	
Dio Implant® Compatible *	AG-UF-RP00	
NeoBiotech® Compatible *	AG-IS-RP00	

Abutment Gauge - 17° Regular Platform

	H	1-5mm
NobelActive™ Compatible *	AG-NA-RP1700	
Osstem® Hiossen® Compatible *	AG-TS-RP1700	
Dentium® Compatible *	AG-DE-RP1700	
MegaGen® AnyOne Compatible *	AG-AO-RP1700	
MegaGen® AnyRidge Compatible *	AG-AR-RP1700	
Astra Tech™TX Compatible *	AG-OS-WP1700	
Dio Implant® Compatible *	AG-UF-RP1700	
NeoBiotech® Compatible *	AG-IS-RP1700	

Abutment Gauge - 30° Regular Platform

	H	1-5mm
NobelActive™ Compatible *	AG-NA-RP3000	
Osstem® Hiossen® Compatible *	AG-TS-RP3000	
Dentium® Compatible *	AG-DE-RP3000	
MegaGen® AnyOne Compatible *	AG-AO-RP3000	
MegaGen® AnyRidge Compatible *	AG-AR-RP3000	
Astra Tech™TX Compatible *	AG-OS-WP3000	
Dio Implant® Compatible *	AG-UF-RP3000	
NeoBiotech® Compatible *	AG-IS-RP3000	

Abutment Gauge Straight Narrow Platform

	H	1-5mm
NobelActive™ Compatible *	AG-NA-NP00	
Osstem® Hiossen® Compatible *	AG-TS-NP00	
Astra Tech™TX Compatible *	AG-OS-SP00	

Abutment Gauge - 17° Narrow Platform

	H	1-5mm
NobelActive™ Compatible *	AG-NA-NP1700	
Osstem® Hiossen® Compatible *	AG-TS-NP1700	
Astra Tech™TX Compatible *	AG-OS-SP1700	

Abutment Gauge - 30° Narrow Platform

	H	1-5mm
NobelActive™ Compatible *	AG-NA-NP3000	
Osstem® Hiossen® Compatible *	AG-TS-NP3000	
Astra Tech™TX Compatible *	AG-OS-SP3000	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.




Internal Hex

Lamina® INTERNAL HEX Connection

 Regular Platform


		3.5mm	4.2mm	5.0mm
8mm		LM-SV-RP3508	LM-SV-RP4208	LM-SV-RP5008
10mm		LM-SV-RP3510	LM-SV-RP4210	LM-SV-RP5010
11.5mm		LM-SV-RP3511	LM-SV-RP4211	LM-SV-RP5011
13mm		LM-SV-RP3513	LM-SV-RP4213	LM-SV-RP5013
15mm		LM-SV-RP3515	LM-SV-RP4215	



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.




Titanium Implant Analog
Regular Platform



Internal Hex
Compatible*

IA-SV-RP


Angled Closed Tray 15°
Regular Platform



Internal Hex
Compatible*

ACTT-SV-RP1500

Angled Closed Tray 25°
Regular Platform




Internal Hex
Compatible*

ACTT-SV-RP2500

Patent Pending US 62/416.368

Clicq Trans™
Angled Closed Tray Impression Coping

Closed Tray Impression Coping - Narrow - Ø 4mm
Regular Platform




Internal Hex
Compatible*

CT-SV-RPSN

CT-SV-RPLN

Included

Closed Tray Impression Coping - Ø 4.5mm
Regular Platform




Internal Hex
Compatible*

CT-SV-RPS

CT-SV-RPL

Open Tray Impression Coping
Regular Platform



Internal Hex
Compatible*

OT-SV-RPN

OT-SV-RP



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Healing Abutment Titanium

Titanium Healing Abutment 3mm Regular Platform

	3mm 3.8mm	3mm 4.5mm	3mm 5.5mm	3mm 6.3mm
H Ø	HA-SV-RP30N	HA-SV-RP30	HA-SV-RP30W	HA-SV-RP30EW
Internal Hex Compatible*				

Titanium Healing Abutment 5mm Regular Platform

	5mm 3.8mm	5mm 4.5mm	5mm 5.5mm	5mm 6.3mm
H Ø	HA-SV-RP50N	HA-SV-RP50	HA-SV-RP50W	HA-SV-RP50EW
Internal Hex Compatible*				

Concave Titanium Healing Abutment 3mm Regular Platform

	3mm 4.5mm
H Ø	HAC-SV-RP30
Internal Hex Compatible*	

Concave Titanium Healing Abutment 5mm Regular Platform

	5mm 4.5mm
H Ø	HAC-SV-RP50
Internal Hex Compatible*	


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex

Temporary Abutment



Temporary Titanium Abutment Regular Platform



PEEK - Straight Abutment Regular Platform

Medical Peek - Classix



PEEK - Anatomic Angled Abutment Regular Platform 15°




PEEK - Anatomic Angled Abutment Regular Platform 25°



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex  **Standard Abutment**
Titanium


 **Straight Abutment - Ø 4.5mm**
Regular Platform



Internal Hex
Compatible*

SA-SV-RP5 SA-SV-RP7 SA-SV-RP9 SA-SV-RP12 SA-SV-RP15

 **Straight Abutment - Wide - Ø 5.5mm**
Regular Platform



Internal Hex
Compatible*

SA-SV-RP9W SA-SV-RP12W


 **Straight Slim Abutment**
Regular Platform



Internal Hex
Compatible*

SSA-SV-RP5 SSA-SV-RP7 SSA-SV-RP9


 **Angled Abutment 15°**
Regular Platform



Internal Hex
Compatible*

TA-SV-RP1500N TA-SV-RP1500 TA-SV-RP1500L


 **Angled Abutment 25°**
Regular Platform



Internal Hex
Compatible*

TA-SV-RP2500 TA-SV-RP2500L

 **Angled Abutment 35°**
Regular Platform



Internal Hex
Compatible*

TA-SV-RP3500

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex  **Straight Abutment**
Titanium


 **Straight Abutment W/Shoulder**
Regular Platform




 **Concave Straight Abutment W/Shoulder**
Regular Platform



 **Anatomic Straight Abutment**
Regular Platform



 **Concave Anatomic Straight Abutment**
Regular Platform



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex




Angled Abutment

Titanium

 Anatomic Angled Abutment - 15°
Regular Platform



 Concave Anatomic Angled Abutment - 15°
Regular Platform





Anatomic Angled Abutment - 25° Regular Platform



 Concave Anatomic Angled Abutment - 25°
Regular Platform




* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Clicq™ Pression Components

Abutment Analog



Medical Peek - Classix

Healing Cap PEEK



Waxing Sleeve



Plastic Transfer



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Transfer Abutment Standard Kit



Contents:

Straight / Angled Abutment
Titanium Screw
Plastic Transfer

Straight Transfer Abutment - Short Regular Platform



Internal Hex Compatible* TAP-SV-RP10S TAP-SV-RP20S TAP-SV-RP30S

Straight Transfer Abutment - Medium Regular Platform



Internal Hex Compatible* TAP-SV-RP10M TAP-SV-RP20M TAP-SV-RP30M

Straight Transfer Abutment - Long Regular Platform




Internal Hex Compatible* TAP-SV-RP10L TAP-SV-RP20L TAP-SV-RP30L

Angled Transfer Abutment - 15° Regular Platform




Internal Hex Compatible* TAP-SV-RP1510 TAP-SV-RP1520 TAP-SV-RP1530

Angled Transfer Abutment - 25° Regular Platform



Internal Hex Compatible* TAP-SV-RP2510 TAP-SV-RP2520 TAP-SV-RP2530

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Professional Kit REF

Internal Hex
Compatible *

LOCW-SV-RPPK



Pop Clicq™ - Components

LOCW-RCV Violet 2.7kg (4 pcs.)	LOCW-RCT Transparent 1.8kg (4 pcs.)	LOCW-RCP Pink 1.2kg (4 pcs.)	LOCW-RCY Yellow 0.6kg (4 pcs.)	LOCW-RCB Black laboratory (4 pcs.)
	Pick-Up Impression Transfer Assembled with Yellow Retentive Cap		Analog for Abutment	Protective Disc
LOCW-MH				OD-TK
LOCW-PU		LOCW-AN	LOCW-PD	<p>Overdenture Tool KIT for Insertion & Extraction of Retentive Caps Suitable for all Overdenture lines</p>

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Pop Clicq™ Mini Kit



Straight Titanium Pop-Clicq™ - Mini Kit Regular Platform



Angled Titanium Pop-Clicq™ - 18° Regular Platform



Angled Titanium Pop-Clicq™ - 30° Regular Platform



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Pop Clicq™ - Bar



Kits REF (Type A)

LOCW-BK2

LOCW-BK2D (Screw Driver Included)

Kits REF (Type B)

LOCW-BK4

LOCW-BK4D (Screw Driver Included)

Pop Clicq™ - Bar Components



Bars *2



(A) Adjustable Metal Housing *2
(B) Adjustable Metal Housing *4
(With a fixed screw)



Pop Clicq TM Bar
(A) Metal Housing Tool \ Adapter *2
(B) Metal Housing Tool \ Adapter *4



(A) Analog for Pop Clicq™ Abutment *2
(B) Analog for Pop Clicq™ Abutment *4



(A) Pick-Up Impression Transfer *2
(B) Pick-Up Impression Transfer *4



Adjusting Screws for
(A) Metal Housing *3 (Incl. spares)
(B) Metal Housing *6 (Incl. spares)

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Ball Attachment Titanium




Ball Attachment Regular Platform



Internal Hex Compatible*

Ball Attachment - Standard Kit Regular Platform



Include:

- BAT-MH (Metal Housing)
- BAT-RCT (Retentive Cap - Transparent)
- LOCN-PD (Protective Disc)



Include:

- BAT-MH (Metal Housing)
- BAT-RCP (Retentive Cap - Pink)
- LOCN-PD (Protective Disc)

Ball Attachment - Components

Retentive Cap

BAT-RCT Transparent Slightly elastic 1.2kg - 1.3kg (4 pcs.)	BAT-RCP Pink Elastic 0.8kg - 0.85kg (4 pcs.)	BAT-RCY Yellow Very elastic 0.5kg - 0.55kg (4 pcs.)	BAT-RCG Green Extremely elastic Up to 0.5kg (4 pcs.)	BAT-RCB Black laboratory (4 pcs.)

Metal Housing



BAT-MH




Protective Disc



LOCN-PD

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Professional Kit REF

Internal Hex
Compatible *

LOCN-SV-RPPK



Q Clicq™ - Components



LOCN-RCV
Violet
2.7kg (4 pcs.)



LOCN-RCT
Transparent
1.8kg (4 pcs.)

Retentive Cap



LOCN-RCP
Pink
1.2kg (4 pcs.)



LOCN-RCY
Yellow
0.6kg (4 pcs.)



LOCN-RCB
Black laboratory (4 pcs.)



Metal Housing
LOCN-MH



Impression Plastic
Transfer

Pick-Up Impression Transfer
Assembled with
Yellow Retentive Cap



LOCN-PU

Analog for
Abutment



LOCN-AN

Titanium Screw



SMUL-SV-1.72

Protective Disc



LOCN-PD



OD-TK

Overture Tool KIT
for Insertion & Extraction
of Retentive Caps Suitable
for all Overture lines

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Q-Clicq™
Mini Kit - Titanium



Straight Titanium Q-Clicq™ Regular Platform



Angled Titanium Q-Clicq™ - 18° Regular Platform



Angled Titanium Q-Clicq™ - 30° Regular Platform



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Professional Kit REF

Internal Hex
Compatible*

MUL-SV-RPPK



Starter Kit REF

Internal Hex
Compatible*

MUL-SV-RPSK

Multi Clicq™ - Components

Abutment Analog



MUL-KS-AN

Closed Tray Plastic Transfer

MUL-KS-PUP
[MUL-KS-TSSV - included]

Castable Plastic Sleeve

MUL-KS-SLVPSC
[MUL-KS-TSSV - included]

Cobalt Sleeve UCLA

MUL-KS-SLVCC
[MUL-KS-TSSV - included]

Open Tray Transfer

MUL-KS-OT
[MUL-KS-SOT - included]

Healing cap

MUL-KS-HA
[MUL-KS-TSSV - included]

Scan Body (Peek)

MUL-KS-SB
[MUL-KS-TSSV - included]

Temporary Titanium Sleeve

MUL-KS-SLVTI
[MUL-KS-TSSV - included]

Rotational Titanium Base

MUL-KS-TBR
[MUL-KS-TSSV - included]

Anti Rotational Titanium Base

MUL-KS-TBT
[MUL-KS-TSSV - included]

Rotational Adhesive Cap

MUL-KS-ADR
[MUL-KS-TSSV - included]

Anti Rotational Adhesive Cap

MUL-KS-AD
[MUL-KS-TSSV - included]

*For Single Unit Restoration.

Prosthetic Titanium Screw



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex



Multi Clicq™ - On4/On6 Mini & Standard Kit



Straight Mini Kit
Regular Platform



H

1mm

2mm

3mm

4mm

5mm

6mm

Internal Hex
Compatible*

MUL-SV-RPKS10K

MUL-SV-RPKS20K

MUL-SV-RPKS30K

MUL-SV-RPKS40K

MUL-SV-RPKS50K

MUL-SV-RPKS60K



Straight Abutment
Regular Platform



1mm

2mm

3mm

4mm

5mm

6mm

H

Internal Hex
Compatible*

MUL-SV-RPKS10

MUL-SV-RPKS20

MUL-SV-RPKS30

MUL-SV-RPKS40

MUL-SV-RPKS50

MUL-SV-RPKS60

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



One Piece Angled 17°
Abutment Mini Kit
Regular Platform



G	H	1mm 2.5mm	2mm 3.5mm
Internal Hex Compatible *		MUL-SV-RPKS1710K	MUL-SV-RPKS1720K



One Piece Angled Abutment - 17°
Standard Kit
Regular Platform

H G	1mm 2.5mm	2mm 3.5mm
Internal Hex Compatible *	MUL-SV-RPKS1710	MUL-SV-RPKS1720

Anti-Rotational - One Piece Angled Abutment - 17°
Standard Kit
Regular Platform

H G	1mm 2.5mm	2mm 3.5mm
Internal Hex Compatible *	MUL-SV-RPKS1710H	MUL-SV-RPKS1720H

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex

Design Patent Pending No. 62154
Patent Pending No. US 62/679,389



Clicq™ - Multi Bar On4/On6 Kit



Clicq™ - Multi Bar On4 Kit Regular Platform



(Screw Driver Included)

Internal Hex
Compatible *

MUL-SV-BK4

MUL-SV-BK4D

Clicq™ - Multi Bar On6 Kit Regular Platform



(Screw Driver Included)

Internal Hex
Compatible *

MUL-SV-BK6

MUL-SV-BK6D

Get a permanent bar for your patient
within minutes



Clicq™ - Multi Bar Components



Bars *2



Pre-Fabricated Permanent
Titanium Adjustable Sleeve *4/6
with Adjusting Micro Screw *8/12



Adjusting Micro Screw *2
(spares)



Prosthetic KS Ti-Screw
On 4 - *6 (Incl. 2 spares)
On 6 - *8 (Incl. 2 spares)




Hand Screw Driver

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Zircocam™ - Straight Anatomic Abutment
Regular Platform




Zircocam™ - Angled Anatomic Abutment 15°
Regular Platform




Zircocam™ - Angled Anatomic Abutment 25°
Regular Platform



Zircotec™ - Straight Anatomic Abutment
Regular Platform




Zircotec™ - Angled Anatomic Abutment 15°
Regular Platform



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Plastic Abutment for Casting
Regular Platform




Internal Hex
Compatible*

PAC-SV-RP

PAC-SV-RPR

Titanium Abutment / UCLA
Regular Platform

(Castable Plastic Sleeve included)



Internal Hex
Compatible*

ABZ-SV-RPTI

ABZ-SV-RPTIR

Chrome Cobalt Abutment / UCLA
Regular Platform

(Castable Plastic Sleeve included)



Internal Hex
Compatible*


ABZ-SV-RPCC

ABZ-SV-RPCCR

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.



Scan Abutment
Regular Platform



Internal Hex
Compatible*


SB-SV-RP

SB-TBC

Sirona CEREC®
Compatible*



Titanium Base
Regular Platform



Internal Hex
Compatible*

TBC-SV-RP

TBC-SV-RPR

Angled Titanium Base 10°
Regular Platform



Internal Hex
Compatible*

TBC-SV-RP1000

TBC-SV-RP1000R

Pre-Milled Titanium Blank
Regular Platform



Internal Hex
Compatible*

PMAM-SV-RP

PMMD-SV-RP

PMNT-SV-RP



Internal Hex
Compatible*

PMAM-SV-RPR

PMMD-SV-RPR

PMNT-SV-RPR

Internal Hex
Compatible*

PMAM-SV-RPR

PMMD-SV-RPR

PMNT-SV-RPR

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Internal Hex  **Abutment Gauge**
Mini Kit




 **Abutment Gauge Kit**
Regular Platform

Internal Hex
Compatible* | **AG-SV-RPK**




 **Abutment Gauge Straight**
Regular Platform




Internal Hex
Compatible* | **AG-SV-RP00**

 **Abutment Gauge - 17°**
Regular Platform



Internal Hex
Compatible* | **AG-SV-RP1700**

 **Abutment Gauge - 30°**
Regular Platform



Internal Hex
Compatible* | **AG-SV-RP3000**



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.





Tools



 Motor Mount for Implant - Self Loading
Regular Platform



	22mm	28mm
NobelActive™ Compatible *	MDIO-NA-RPS	MDIO-NA-RPL
Osstem® Hiossen® Compatible *	MDIO-TS-RPS	MDIO-TS-RPL
Dentium® Compatible *	MDIO-DE-RPS	MDIO-DE-RPL
MegaGen® AnyOne Compatible *	MDIO-AO-RPS	MDIO-AO-RPL
MegaGen® AnyRidge Compatible *	MDIO-AR-RPS	MDIO-AR-RPL
Astra Tech™TX Compatible *	MDIO-OS-WPS	MDIO-OS-WPL
Dio Implant® Compatible *	MDIO-UF-RPS	MDIO-UF-RPL
NeoBiotech® Compatible *	MDIO-IS-RPS	MDIO-IS-RPL
Internal Hex Compatible *	MDIO-SV-RPS	MDIO-SV-RPL


 Motor Mount for Implant - Self Loading
Narrow Platform



	22mm	28mm
NobelActive™ Compatible *	MDIO-NA-NPS	MDIO-NA-NPL
Osstem® Hiossen® Compatible *	MDIO-TS-NPS	MDIO-TS-NPL
Astra Tech™TX Compatible *	MDIO-OS-SPS	MDIO-OS-SPL


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Ratchet Driver for Implant - Self Loading
Regular Platform




NobelActive™ Compatible*	RDIO-NA-RPXS	RDIO-NA-RPS	RDIO-NA-RPL
Osstem® Hiossen® Compatible*	RDIO-TS-RPXS	RDIO-TS-RPS	RDIO-TS-RL
Dentium® Compatible*	RDIO-DE-RPXS	RDIO-DE-RPS	RDIO-DE-RPL
MegaGen® AnyOne Compatible*	RDIO-AO-RPXS	RDIO-AO-RPS	RDIO-AO-RPL
MegaGen® AnyRidge Compatible*	RDIO-AR-RPXS	RDIO-AR-RPS	RDIO-AR-RPL
Astra Tech™TX Compatible*	RDIO-OS-WPXS	RDIO-OS-WPS	RDIO-OS-WPL
Dio Implant® Compatible*	RDIO-UF-RPXS	RDIO-UF-RPS	RDIO-UF-RPL
NeoBiotech® Compatible*	RDIO-IS-RPXS	RDIO-IS-RPS	RDIO-IS-RPL
Internal Hex Compatible*	RDIO-SV-RPXS	RDIO-SV-RPS	RDIO-SV-RPL

Ratchet Driver for Implant - Self Loading
Regular Platform




NobelActive™ Compatible*	TDIO-NA-RPS	TDIO-NA-RPL
Osstem® Hiossen® Compatible*	TDIO-TS-RPS	TDIO-TS-RPL
Dentium® Compatible*	TDIO-DE-RPS	TDIO-DE-RPL
MegaGen® AnyOne Compatible*	TDIO-AO-RPS	TDIO-AO-RPL
MegaGen® AnyRidge Compatible*	TDIO-AR-RPS	TDIO-AR-RPL
Astra Tech™TX Compatible*	TDIO-OS-WPS	TDIO-OS-WPL
Dio Implant® Compatible*	TDIO-UF-RPS	TDIO-UF-RPL
NeoBiotech® Compatible*	TDIO-IS-RPS	TDIO-IS-RPL
Internal Hex Compatible*	TDIO-SV-RPS	TDIO-SV-RPL

Ratchet Driver for Implant - Self Loading
Narrow Platform



NobelActive™ Compatible*	RDIO-NA-NPXS	RDIO-NA-NPS	RDIO-NA-NPL
Osstem® Hiossen® Compatible*	RDIO-TS-NPXS	RDIO-TS-NPS	RDIO-TS-NPL
Astra Tech™TX Compatible*	RDIO-OS-SPXS	RDIO-OS-SPS	RDIO-OS-SPL

Ratchet Driver for Implant - Self Loading
Narrow Platform



NobelActive™ Compatible*	TDIO-NA-NPS	TDIO-NA-NPL
Osstem® Hiossen® Compatible*	TDIO-TS-NPS	TDIO-TS-NPL
Astra Tech™TX Compatible*	TDIO-OS-SPS	TDIO-OS-SPL

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Tools For Screw




Motor Mount for Screw Regular Platform		
	L	
	22mm	
	28mm	
NobelActive™ Compatible*	MDS-NA-S	MDS-NA-L
Osstem® Hiossen® Compatible*	MDS-TS-S	MDS-TS-L
Dentium® Compatible*	MDS-DE-S	MDS-DE-L
MegaGen® AnyOne Compatible*	MDS-AO-S	MDS-AO-L
MegaGen® AnyRidge Compatible*	MDS-AR-S	MDS-AR-L
Astra Tech™TX Compatible*	MDS-OS-S	MDS-OS-L
Dio Implant® Compatible*	MDS-UF-S	MDS-UF-L
NeoBiotech® Compatible*	MDS-IS-S	MDS-IS-L
Internal Hex Compatible*	MDS-SV-S	MDS-SV-L

Ratchet Driver for Screw Regular Platform		
	4	
	L	
	12mm	
	18mm	
NobelActive™ Compatible*	TDS-NA-RPS	TDS-NA-RPL
Osstem® Hiossen® Compatible*	TDS-TS-RPS	TDS-TS-RPL
Dentium® Compatible*	TDS-DE-RPS	TDS-DE-RPL
MegaGen® AnyOne Compatible*	TDS-AO-RPS	TDS-AO-RPL
MegaGen® AnyRidge Compatible*	TDS-AR-RPS	TDS-AR-RPL
Astra Tech™TX Compatible*	TDS-OS-WPS	TDS-OS-WPL
Dio Implant® Compatible*	TDS-UF-RPS	TDS-UF-RPL
NeoBiotech® Compatible*	TDS-IS-RPS	TDS-IS-RPL
Internal Hex Compatible*	TDS-SV-RPS	TDS-SV-RPL

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Tools For Screw



Ratchet Driver for Screw
Regular Platform



NobelActive™ Compatible *	RDS-NA-S	RDS-NA-L		Ø 8mm
Osstem® Hiossen® Compatible *	RDS-TS-S	RDS-TS-L		
Dentium® Compatible *	RDS-DE-S	RDS-DE-L		
MegaGen® AnyOne Compatible *	RDS-AO-S	RDS-AO-L		
MegaGen® AnyRidge Compatible *	RDS-AR-S	RDS-AR-L		Ø 7mm
Astra Tech™TX Compatible *	RDS-OS-S	RDS-OS-L		
Dio Implant® Compatible *	RDS-UF-S	RDS-UF-L		
NeoBiotech® Compatible *	RDS-IS-S	RDS-IS-L		
Internal Hex Compatible *	MDS-SV-S	MDS-SV-L		

 Hand Screw Driver
Regular Platform




NobelActive™ Compatible *	HSDA-NA-S	HSDA-NA-L	HSDA-NA-XL
Osstem® Hiossen® Compatible *	HSDA-TS-S	HSDA-TS-L	HSDA-TS-XL
Dentium® Compatible *	HSDA-DE-S	HSDA-DE-L	HSDA-DE-XL
MegaGen® AnyOne Compatible *	HSDA-AO-S	HSDA-AO-L	HSDA-AO-XL
MegaGen® AnyRidge Compatible *	HSDA-AR-S	HSDA-AR-L	HSDA-AR-XL
Astra Tech™TX Compatible *	HSDA-OS-S	HSDA-OS-L	HSDA-OS-XL
Dio Implant® Compatible *	HSDA-UF-S	HSDA-UF-L	HSDA-UF-XL
NeoBiotech® Compatible *	HSDA-IS-S	HSDA-IS-L	HSDA-IS-XL
Internal Hex Compatible *	HSDA-SV-S	HSDA-SV-L	HSDA-SV-XL


* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Tools
Drills


 **Drill Extender**
External Irrigation



 **Lance Drill**
External Irrigation



 **Dual Drill Ø 2.0mm**
External Irrigation



 **Step drills - Ø 2.0mm**
External Irrigation



 **Step drills - Ø 2.7mm**
External Irrigation



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Tools
Drills

 Step drills - Ø 3.1mm
External Irrigation

		8mm	10mm	11.5mm	13mm	15mm
L		SD-3.1-08	SD-3.1-10	SD-3.1-11	SD-3.1-13	SD-3.1-15
REF						

 Step drills - Ø 3.7mm
External Irrigation

		8mm	10mm	11.5mm	13mm	15mm
H (mm)		SD-3.7-8	SD-3.7-10	SD-3.7-11	SD-3.7-13	SD-3.7-15
REF						


 Step drills - Ø 4.4mm
External Irrigation

		8mm	10mm	11.5mm	13mm
L		SD-4.4-8	SD-4.4-10	SD-4.4-11	SD-4.4-13
REF					

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.


Tools
Drills

 **Twist Step Drills**
External Irrigation


					
Ø1	2.4mm	2.7mm	3.1mm	3.7mm	4.4mm
Ø2	2.0mm	2.4mm	2.7mm	3.1mm	3.7mm
REF	TSD-2.4-15	TSD-2.7-15	TSD-3.1-15	TSD-3.7-15	TSD-4.4-15

 **Stoppers for Standard Twist Drills**


				
REF	DSTD-13	DSTD-11.5	DSTD-10	DSTD-8



 **Hard Bone Drills**
Regular Platform

			
Ø1	3.5mm	5.0mm	
Ø2	3.1mm	4.2mm	
REF	HBD-3.5	HBD-5.0	

 **Trephine Burs**
Regular Platform


			
Ø1	3mm	4mm	
Ø2	4mm	5mm	
REF	TB-3.0-4.0	TB-4.0-5.0	

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Tools



RW-UH Ratchet Wrench
Black Marking w/6.35mm Hex



RW-UTS Ratchet Wrench
Black Marking w/Square TS 4mm




HDR Hand Driver suitable for
Squar Ratchet Driver
(Ratchet driver not Incl.)



RT-45H Torque Ratchet 0-45 Ncm
Black Marking w/6.35 Hex



RT-45TS Torque Ratchet 0-45 Ncm
Black Marking w/Square TS 4mm



PPN
Parallel Pin Guide



RT-70ST Torque Ratchet 0-70 Ncm
Black Marking, Ø - 7mm



RT-70TSSQ Torque Ratchet 0-70 Ncm
Black Marking, Square TS 4mm



MM-ADP-8
Motor Mount Adapter
w/Ball Friction 8mm

MM-ADP-7
Motor Mount Adapter
w/Ball Friction 7mm



IPR-S
Implant Deep Probe - Single End



TS-ADP-8
Ratchet Adapter
from Square 4mm
to Round Ø 8mm

TS-ADP-7
Ratchet Adapter
from Square 4mm
to Round Ø 7mm



IPR-D
Implant Deep Probe

* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Surgical Tool Box



REF RSB-AR-SK

Surgical Tool Box - Artic Components

Lance Drill:
LD-2.0-06

Torque Ratchet 0-70 Ncm:
RT-70ST


Drill Extender:
DEXT

Tools for Implants:

Twist Step Drills:
TSD-2.4-15
TSD-2.7-15
TSD-3.1-15
TSD-3.7-15
TSD-4.4-15

Ratchet Driver for Implant:

RDIO-NA-RPS
RDIO-NA-RPL
RDIO-NA-NPS
RDIO-NA-NPL



Hard Bone Drills:
HBD-3.5
HBD-5.0

Motor Mount for Implant:

MDIO-NA-RPS
MDIO-NA-RPL
MDIO-NA-NPS
MDIO-NA-NPL

Stoppers for Twist Drills:
DSTD-13 *2
DSTD-11.5 *2
DSTD-10 *2
DSTD-8 *2

Tools for Screw:

Ratchet Driver for Screw:
RDS-SV-S
RDS-SV-L

Parallel Pin:
PPN *3

Motor Mount for Screw:
MDS-SV-L

Hand Driver for Screw:
HSDA-SV-S
HSDA-SV-L



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Surgical Tool Box



REF RSB-FR-SK

Surgical Tool Box - ios3/Fortis Components

Lance Drill:
LD-2.0-06

Torque Ratchet 0-70 Ncm:
RT-70ST


Drill Extender:
DEXT

Tools for Implants:

Twist Step Drills:
TSD-2.4-15
TSD-2.7-15
TSD-3.1-15
TSD-3.7-15
TSD-4.4-15

Ratchet Driver for Implant:

RDIO-TS-RPS
RDIO-TS-RPL
RDIO-OS-SPS
RDIO-OS-SPL



Hard Bone Drills:
HBD-3.5
HBD-5.0

Motor Mount for Implant:

MDIO-TS-RPS
MDIO-TS-RPL
MDIO-OS-SPS
MDIO-OS-SPL

Stoppers for Twist Drills:
DSTD-13 *2
DSTD-11.5 *2
DSTD-10 *2
DSTD-8 *2

Tools for Screw:

Ratchet Driver for Screw:
RDS-SV-S
RDS-SV-L

Parallel Pin:
PPN *3

Motor Mount for Screw:
MDS-SV-L

Hand Driver for Screw:
HSDA-SV-S
HSDA-SV-L



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.

Surgical Tool Box



REF RSB-LM-SK

Surgical Tool Box - Lamina Components

Lance Drill:
LD-2.0-06

Torque Ratchet 0-70 Ncm:
RT-70ST

Drill Extender:
DEXT

Tools for Implants:
Ratchet Driver for Implant:

Twist Step Drills:
TSD-2.4-15
TSD-2.7-15
TSD-3.1-15
TSD-3.7-15
TSD-4.4-15

RDIO-SV-RPS
RDIO-SV-RPL

Hard Bone Drills:
HBD-3.5
HBD-5.0

Tools for Screw:


Stoppers for Twist Drills:
DSTD-13 *2
DSTD-11.5 *2
DSTD-10 *2
DSTD-8 *2

Ratchet Driver for Screw:
RDS-SV-S
RDS-SV-L

Parallel Pin:
PPN *3

Motor Mount for Screw:
MDS-SV-L

Hand Driver for Screw:
HSDA-SV-S
HSDA-SV-L



* Registered trademarks are owned by their respective companies and there is no commercial relationship with Edison Medical. Its mention is made for proper identification purpose of the dental implants with which the Edison Medical abutments can be used.